



**HEALTH MISSION
INSTRUMENTS**

Catalog

General Surgical Instruments

2012-2013







Catalog

Foreword

Procedure	Page No
General	1
Plastic and Maxillo Facial	78
Ear, Nose and Throat	110
Thoracic	170
Orthopaedic	196
Gynaecology and Obstetrics	222
Neurology	236
Proctology	257
Urology	266
Ophthalmic	279


General Artery Forceps

11000-18 Adson (Fraser) Artery Forcep, curved, 1/2 teeth, 180mm

11010-18 Adson (Fraser) Artery Forcep, curved, 180mm

11020-18 Adson (Fraser) Artery Forcep, straight, 1/2 teeth, 180mm

11030-18 Adson (Fraser) Artery Forcep, straight, 180mm



11040-35 Dieffenbach Artery Forcep, curved, 35mm

11040-50 Dieffenbach Artery Forcep, curved, 50mm

11050-35 Dieffenbach Artery Forcep, straight, 35mm


11050-50 Dieffenbach Artery Forcep, straight, 50mm




11060-12 Dunhill Artery Forcep, curved, box joint, 125mm

11060-18 Dunhill Artery Forcep, curved, box joint, 180mm

11070-12 Dunhill Artery Forcep, straight, box joint, 125mm




11080-18 Fraser Kelly Artery Forcep, curved, 180mm




11090-14 Grey Turner Artery Forcep, curved, 145mm


11090-18 Grey Turner Artery Forcep, curved, 180mm




General Artery Forceps




- 11100-14 Grey Turner Artery Forcep, straight, 145mm
11100-18 Grey Turner Artery Forcep, straight, 180mm




- 11110-12 Halstead Mosquito Artery Forcep, curved, 1/2 teeth, 125mm
11130-12 Halstead Mosquito Artery Forcep, straight, 1/2 teeth, 125mm



- 11120-09 Halstead Mosquito Artery Forcep, curved, serrated jaws, 95mm
11120-12 Halstead Mosquito Artery Forcep, curved, serrated jaws, 125mm
11140-09 Halstead Mosquito Artery Forcep, straight, serrated jaws, 95mm
11140-12 Halstead Mosquito Artery Forcep, straight, serrated jaws, 125mm



- 11150-20 Heiss Artery Forcep, curved, 200mm
11160-20 Heiss Artery Forcep, slightly curved, 200mm
11170-20 Heiss Artery Forcep, straight, 200mm




- 11180-14 Kelly-Crile Artery Forcep, curved, half serrated jaws, 145mm
11180-16 Kelly-Crile Artery Forcep, curved, half serrated jaws, 160mm
11190-14 Kelly-Crile Artery Forcep, straight, half serrated jaws, 145mm
11190-16 Kelly-Crile Artery Forcep, straight, half serrated jaws, 160mm




General Artery Forceps


- 11200-14 Kilner Artery Forcep, curved, serrated jaws, 140mm
11210-14 Kilner Artery Forcep, straight, serrated jaws, 140mm




- 11220-12 Kocher Artery Forcep, curved, 1/2 teeth, 125mm
11220-14 Kocher Artery Forcep, curved, 1/2 teeth, 140mm
11220-15 Kocher Artery Forcep, curved, 1/2 teeth, 150mm
11220-18 Kocher Artery Forcep, curved, 1/2 teeth, 180mm
11220-20 Kocher Artery Forcep, curved, 1/2 teeth, 200mm
11220-23 Kocher Artery Forcep, curved, 1/2 teeth, 230mm




- 11230-12 Kocher Artery Forcep, straight, 1/2 teeth, 125mm
11230-14 Kocher Artery Forcep, straight, 1/2 teeth, 140mm
11230-15 Kocher Artery Forcep, straight, 1/2 teeth, 150mm
11230-18 Kocher Artery Forcep, straight, 1/2 teeth, 180mm
11230-20 Kocher Artery Forcep, straight, 1/2 teeth, 200mm
11230-23 Kocher Artery Forcep, straight, 1/2 teeth, 230mm




- 11240-18 Mayo Oschner Artery Forcep, curved, 185mm
11250-18 Mayo Oschner Artery Forcep, straight, 185mm



- 11260-15 Mayo Artery Forcep, curved, 150mm




General Artery Forceps




11270-14 Mixer Baby Artery Forcep, curved, half serrated jaws, 140mm

11270-18 Mixer Baby Artery Forcep, curved, half serrated jaws, 180mm




11275-15 Mixer Artery Forcep, curved, serrated jaws, 150mm

11275-18 Mixer Artery Forcep, curved, serrated jaws, 180mm




11275-22 Mixer Artery Forcep, curved, half serrated jaws, 220mm



11280-14 Moynihan Artery Forcep, curved, half serrated jaws, 145mm

11280-18 Moynihan Artery Forcep, curved, half serrated jaws, 180mm




11290-14 Moynihan Artery Forcep, straight, half serrated jaws, 145mm

11290-18 Moynihan Artery Forcep, straight, half serrated jaws, 180mm




General Artery Forceps


- 11300-15 O'Shaughnessy Artery Forcep, curved, 150mm
 11300-18 O'Shaughnessy Artery Forcep, curved, 180mm
 11300-20 O'Shaughnessy Artery Forcep, curved, 200mm




- 11310-14 Oschner Artery Forcep, curved, 1/2 teeth, 140mm
 11310-16 Oschner Artery Forcep, curved, 1/2 teeth, 165mm
 11310-18 Oschner Artery Forcep, curved, 1/2 teeth, 180mm
 11310-20 Oschner Artery Forcep, curved, 1/2 teeth, 200mm
 11310-22 Oschner Artery Forcep, curved, 1/2 teeth, 220mm
 11310-24 Oschner Artery Forcep, curved, 1/2 teeth, 240mm




- 11320-14 Oschner Artery Forcep, straight, 1/2 teeth, 140mm
 11320-16 Oschner Artery Forcep, straight, 1/2 teeth, 165mm
 11320-18 Oschner Artery Forcep, straight, 1/2 teeth, 180mm
 11320-20 Oschner Artery Forcep, straight, 1/2 teeth, 200mm
 11320-22 Oschner Artery Forcep, straight, 1/2 teeth, 220mm
 11320-24 Oschner Artery Forcep, straight, 1/2 teeth, 240mm




- 11330-14 Rochester-Pean Artery Forcep, curved, 140mm
 11330-16 Rochester-Pean Artery Forcep, curved, 160mm
 11330-18 Rochester-Pean Artery Forcep, curved, 180mm
 11330-20 Rochester-Pean Artery Forcep, curved, 200mm
 11330-22 Rochester-Pean Artery Forcep, curved, 220mm
 11330-24 Rochester-Pean Artery Forcep, curved, 240mm
 11330-26 Rochester-Pean Artery Forcep, curved, 260mm
 11340-14 Rochester-Pean Artery Forcep, straight, 140mm
 11340-16 Rochester-Pean Artery Forcep, straight, 160mm
 11340-18 Rochester-Pean Artery Forcep, straight, 180mm
 11340-20 Rochester-Pean Artery Forcep, straight, 200mm
 11340-22 Rochester-Pean Artery Forcep, straight, 220mm
 11340-24 Rochester-Pean Artery Forcep, straight, 240mm
 11340-26 Rochester-Pean Artery Forcep, straight, 260mm




General Artery Forceps




11345-10 Ronald Edwards Artery Forcep, curved, screw joint, shape A, 215mm



11345-20 Ronald Edwards Artery Forcep, curved, screw joint, shape B, 215mm



11345-30 Ronald Edwards Artery Forcep, curved, screw joint, shape C, 215mm



11350-12 Spencer-Wells Artery Forcep, curved, serrated jaws, 125mm

11350-15 Spencer-Wells Artery Forcep, curved, serrated jaws, 150mm

11350-18 Spencer-Wells Artery Forcep, curved, serrated jaws, 180mm


11350-20 Spencer-Wells Artery Forcep, curved, serrated jaws, 200mm

11350-23 Spencer-Wells Artery Forcep, curved, serrated jaws, 230mm



General Artery Forceps

- 11360-12 Spencer-Wells Artery Forcep, straight, serrated jaws, 125mm
- 11360-15 Spencer-Wells Artery Forcep, straight, serrated jaws, 150mm
- 11360-18 Spencer-Wells Artery Forcep, straight, serrated jaws, 180mm
- 11360-20 Spencer-Wells Artery Forcep, straight, serrated jaws, 200mm
- 11360-23 Spencer-Wells Artery Forcep, straight, serrated jaws, 230mm




Brushes

- 11800-28 Wire Brush, for cleaning instruments, 280mm




Cholecystectomy Forceps


- 11900-23 Geary-Grant Cholecystectomy Forcep, large, non-toothed, longitudinal serrations, 230mm
- 11910-23 Geary-Grant Cholecystectomy Forcep, small, non-toothed, longitudinal serrations, 230mm



- 11920-19 Lahey Cholecystectomy Forcep, curved, longitudinal serrations, 195mm
- 11920-23 Lahey Cholecystectomy Forcep, curved, longitudinal serrations, 230mm




- 11930-21 Moynihan Cholecystectomy Forcep, curved on flat, box joint, 215mm



General Clamps



- 12010-10 De Martel Anastomosis Clamp, closing device only
 12010-20 De Martel Anastomosis Clamp, segment clamp only
 12010-30 De Martel Anastomosis Clamp, segment only
 12010-40 De Martel Anastomosis Clamp, tommy bar only
 12010-90 De Martel Anastomosis Clamp, complete with 3 segments, tommy bar and segment clamp



- 12020-17 Doyen Intestinal Clamp, curved, 170mm
 12020-23 Doyen Intestinal Clamp, curved, 230mm




- 12030-17 Doyen Intestinal Clamp, straight, 170mm
 12030-23 Doyen Intestinal Clamp, straight, 230mm



- 12040-23 Guard For Parker Kerr Clamp, fenestrated on inner curve, to suit 230mm clamp
 12040-25 Guard For Parker Kerr Clamp, fenestrated on inner curve, to suit 255mm clamp



- 12050-23 Guard For Parker Kerr Clamp, fenestrated on outer curve, to suit 230mm clamp
 12050-25 Guard For Parker Kerr Clamp, fenestrated on outer curve, to suit 255mm clamp




- 12060-10 Kocher Intestinal Clamp, curved, longitudinal serrations, 225mm, 100mm blades
 12060-12 Kocher Intestinal Clamp, curved, longitudinal serrations, 250mm, 125mm blades
 12060-15 Kocher Intestinal Clamp, curved, longitudinal serrations, 275mm, 150mm blades
 12070-10 Kocher Intestinal Clamp, straight, longitudinal serrations, 225mm, 100mm blades
 12070-12 Kocher Intestinal Clamp, straight, longitudinal serrations, 250mm, 125mm blades
 12070-15 Kocher Intestinal Clamp, straight, longitudinal serrations, 275mm, 150mm blades




General Clamps

- 12080-21 Lane Intestinal Clamp, straight, longitudinal serrations, 215mm, 75mm blades




- 12090-30 Lane Twin Clamp, curved, longitudinal serrations, large, 300mm

- 12100-30 Lane Twin Clamp, straight, longitudinal serrations, large, 300mm



- 12110-20 Lang Stevenson Clamp, keyed blades, 200mm, 80mm blades




- 12120-23 Parker Kerr Clamp, curved, screw joint, 230mm

- 12120-25 Parker Kerr Clamp, curved, screw joint, 255mm


- 12130-23 Parker Kerr Clamp, straight, screw joint, 230mm

- 12130-25 Parker Kerr Clamp, straight, screw joint, 255mm




- 12140-20 Payr Intestinal Clamp, lever action, longitudinal serrations, 200mm, 60mm blades

- 12140-25 Payr Intestinal Clamp, lever action, longitudinal serrations, 255mm, 100mm blades




General Clamps




12150-17 Schoemaker-Ogilvie Clamp, straight, hollow blades, 170mm

12150-21 Schoemaker-Ogilvie Clamp, straight, hollow blades, 215mm



12160-14 Tubing Clamp, 145mm

12160-18 Tubing Clamp, 180mm




12170-10 Zachary Cope Anastomosis Clamp, closing device only

12170-20 Zachary Cope Anastomosis Clamp, segment only

12170-30 Zachary Cope Anastomosis Clamp, tommy bar only

12170-90 Zachary Cope Anastomosis Clamp, complete with 3 segments and tommy bar

Clamps - Atraumatic




12200-18 Doyen Atraumatic Clamp, curved, 180mm




General Clamps-Atraumatic


12200-23 Doyen Atraumatic Clamp, curved, 230mm



12210-18 Doyen Atraumatic Clamp, straight, 180mm



12210-23 Doyen Atraumatic Clamp, straight, 230mm




Clip Forceps

12300-14 Comyns Berkeley Clip Forcep, applying, serrated jaws, 140mm




12320-12 Michel Clip Forcep, applying and extracting combined, 125mm



12340-12 Michel Clip Forcep, applying, 125mm




General Clip Forceps



12360-12 Michel Clip Forcep, clip extracting, 120mm

Decompressors



12400-50 Savage Intestinal Decompressor, s/s, trocar and cannula, 8mm o.d., 500mm

Depressors




12500-25 Sargent Intestinal Depressor, 250mm, 65mm wide to 40mm at narrow end



12520-45 St. Barts Intestinal Depressor, 45mm x 28mm

Diathermy Instruments




13060-20 Diathermy Quiver, plastic with ring, 200mm, 65mm diameter



General Dilators


- 13200-03 Bake Bile Duct Dilator, s/s, malleable stems, 200mm, 3mm diameter head
- 13200-04 Bake Bile Duct Dilator, s/s, malleable stems, 200mm, 4mm diameter head
- 13200-05 Bake Bile Duct Dilator, s/s, malleable stems, 200mm, 5mm diameter head
- 13200-06 Bake Bile Duct Dilator, s/s, malleable stems, 200mm, 6mm diameter head
- 13200-07 Bake Bile Duct Dilator, s/s, malleable stems, 200mm, 7mm diameter head
- 13200-08 Bake Bile Duct Dilator, s/s, malleable stems, 200mm, 8mm diameter head
- 13200-09 Bake Bile Duct Dilator, s/s, malleable stems, 200mm, 9mm diameter head
- 13200-10 Bake Bile Duct Dilator, s/s, malleable stems, 200mm, 10mm diameter head
- 13200-11 Bake Bile Duct Dilator, s/s, malleable stems, 200mm, 11mm diameter head
- 13200-90 Bake Bile Duct Dilator, s/s, malleable stems, 200mm, set of 9



- 13205-03 Clarke Bile Duct Dilator, s/s, double ended, spiral wire stem, 230mm, 3/4mm
- 13205-05 Clarke Bile Duct Dilator, s/s, double ended, spiral wire stem, 230mm, 5/6mm
- 13205-07 Clarke Bile Duct Dilator, s/s, double ended, spiral wire stem, 230mm, 7/8mm



- 13220-04 Maingot Bile Duct Dilator, s/s, malleable stems, double ended, 180mm, 3mm and 4mm
- 13220-06 Maingot Bile Duct Dilator, s/s, malleable stems, double ended, 180mm, 5mm and 6mm
- 13220-08 Maingot Bile Duct Dilator, s/s, malleable stems, double ended, 180mm, 7mm and 8mm
- 13220-90 Maingot Bile Duct Dilator, s/s, malleable stems, double ended, 180mm, set of 3



Directors

- 13300-14 Brodie Director, probe-ended, s/p, 140mm
- 13300-16 Brodie Director, probe-ended, s/p, 165mm
- 13300-20 Brodie Director, probe-ended, s/p, 200mm
- 13310-14 Brodie Director, probe-ended, silver, 140mm
- 13310-16 Brodie Director, probe-ended, silver, 165mm
- 13310-20 Brodie Director, probe-ended, silver, 200mm
- 13320-14 Brodie Director, probe-ended, s/s, 140mm
- 13320-16 Brodie Director, probe-ended, s/s, 165mm
- 13320-20 Brodie Director, probe-ended, s/s, 200mm



- 13340-15 Key Hernia Director, 155mm



General Dissecting Forceps



- 13400-12 Adson Dissecting Forcep, 1/2 teeth, 125mm
 13410-12 Adson Dissecting Forcep, 1/2 teeth, fine, 125mm
 13410-18 Adson Dissecting Forcep, 1/2 teeth, fine, 180mm
 13415-18 Adson Dissecting Forcep, 1/2 teeth, ultra light, 180mm



- 13420-12 Adson Dissecting Forcep, serrated, 125mm
 13430-12 Adson Dissecting Forcep, serrated, fine, 125mm
 13430-18 Adson Dissecting Forcep, serrated, fine, 180mm
 13435-18 Adson Dissecting Forcep, serrated, ultra light, 180mm



- 13440-12 Block End Dissecting Forcep, round serrated jaws, narrow flutes, 125mm
 13440-15 Block End Dissecting Forcep, round serrated jaws, narrow flutes, 150mm
 13440-18 Block End Dissecting Forcep, round serrated jaws, narrow flutes, 180mm
 13440-20 Block End Dissecting Forcep, round serrated jaws, narrow flutes, 200mm
 13440-25 Block End Dissecting Forcep, round serrated jaws, narrow flutes, 250mm



- 13441-15 Block End Dissecting Forcep, fine serrated jaws, narrow flutes, 150mm
 13441-18 Block End Dissecting Forcep, fine serrated jaws, narrow flutes, 180mm
 13441-20 Block End Dissecting Forcep, fine serrated jaws, narrow flutes, 200mm



- 13442-12 Block End Dissecting Forcep, serrated jaws, broad flutes, 125mm
 13442-15 Block End Dissecting Forcep, serrated jaws, broad flutes, 150mm
 13442-18 Block End Dissecting Forcep, serrated jaws, broad flutes, 180mm
 13442-20 Block End Dissecting Forcep, serrated jaws, broad flutes, 200mm
 13442-25 Block End Dissecting Forcep, serrated jaws, broad flutes, 250mm
 13442-30 Block End Dissecting Forcep, serrated jaws, broad flutes, 300mm



- 13445-23 Bickford Dissecting Forcep, 1/2 teeth, 230mm



General Dissecting Forceps

13446-23 Bickford Dissecting Forcep, serrated jaws, 230mm




13450-18 Bonney Dissecting Forcep, 1/2 teeth, 180mm

13460-18 Bonney Dissecting Forcep, serrated, 180mm




13465-14 Braithwaite Dissecting Forcep, 1/2 teeth, 140mm




13470-14 Canadian Dissecting Forcep, 2/3 teeth, 140mm

13470-15 Canadian Dissecting Forcep, 2/3 teeth, 150mm



13480-14 Canadian Dissecting Forcep, serrated jaws, 140mm

13480-15 Canadian Dissecting Forcep, serrated jaws, 150mm



13490-10 Continental Pattern Dissecting Forcep, 1/2 teeth, 105mm

13490-12 Continental Pattern Dissecting Forcep, 1/2 teeth, 125mm

13490-15 Continental Pattern Dissecting Forcep, 1/2 teeth, 150mm

13490-16 Continental Pattern Dissecting Forcep, 1/2 teeth, 160mm

13490-17 Continental Pattern Dissecting Forcep, 1/2 teeth, 175mm

13490-20 Continental Pattern Dissecting Forcep, 1/2 teeth, 200mm

13490-25 Continental Pattern Dissecting Forcep, 1/2 teeth, 250mm

13490-30 Continental Pattern Dissecting Forcep, 1/2 teeth, 300mm



13510-20 Emmett Dissecting Forcep, serrated jaws, 200mm




General Dissecting Forceps



13520-15 Gillies Dissecting Forcep, 1/2 teeth, 150mm



13521-15 Gillies Dissecting Forcep, serrated jaws, 150mm



13530-14 Jeans Dissecting Forcep, 1/2 teeth, 140mm

13530-18 Jeans Dissecting Forcep, 1/2 teeth, 180mm

13540-18 Jeans Dissecting Forcep, 2/3 teeth, 180mm

13550-14 Jeans Dissecting Forcep, serrated jaws, 140mm



13560-12 Lanes Dissecting Forcep, 1/2 teeth, 125mm

13560-14 Lanes Dissecting Forcep, 1/2 teeth, 140mm

13560-18 Lanes Dissecting Forcep, 1/2 teeth, 180mm

13560-20 Lanes Dissecting Forcep, 1/2 teeth, 200mm



13570-12 Lanes Dissecting Forcep, 2/3 teeth, 125mm

13570-14 Lanes Dissecting Forcep, 2/3 teeth, 140mm

13570-18 Lanes Dissecting Forcep, 2/3 teeth, 180mm

13570-20 Lanes Dissecting Forcep, 2/3 teeth, 200mm



13580-14 Leedham Green Dissecting Forcep, 2/3 teeth, 145mm



13590-14 Leedham Green Dissecting Forcep, 3/4 teeth, 145mm



13600-28 Maingot Dissecting Forcep, straight, serrated tips, fenestrated grip, 280mm



General Dissecting Forceps

13610-15 McIndoe Dissecting Forcep, 1/2 teeth, 150mm



13620-15 McIndoe Dissecting Forcep, serrated, 150mm



13640-18 Ramsey Dissecting Forcep, 1/2 teeth, 180mm

13650-18 Ramsey Dissecting Forcep, 2/3 teeth, 180mm



13660-18 Ramsey Dissecting Forcep, serrated, 180mm



13670-29 St.Marks Dissecting Forcep, 1/2 teeth, 290mm

13680-29 St.Marks Dissecting Forcep, serrated, 290mm




13690-12 Treves Dissecting Forcep, 1/2 teeth, 125mm

13690-15 Treves Dissecting Forcep, 1/2 teeth, 150mm

13690-18 Treves Dissecting Forcep, 1/2 teeth, 180mm

13690-20 Treves Dissecting Forcep, 1/2 teeth, 200mm

13690-25 Treves Dissecting Forcep, 1/2 teeth, 250mm



13700-12 Turn-Over-End Dissecting Forcep, serrated jaws, broad flutes, 125mm

13700-15 Turn-Over-End Dissecting Forcep, serrated jaws, broad flutes, 150mm



13710-15 Waugh Dissecting Forcep, 1/2 teeth, 150mm

13710-18 Waugh Dissecting Forcep, 1/2 teeth, 180mm

13710-20 Waugh Dissecting Forcep, 1/2 teeth, 200mm

13720-15 Waugh Dissecting Forcep, serrated jaws, 150mm

13720-18 Waugh Dissecting Forcep, serrated jaws, 180mm

13720-20 Waugh Dissecting Forcep, serrated jaws, 200mm



General Dissectors



13900-19 Durham Dissector, double ended, 190mm



13910-16 Kocher Thyroid Dissector, curved, 160mm



13920-16 Kocher Thyroid Dissector, curved, 165mm



13930-18 Kocher Thyroid Dissector, curved, 180mm



13940-19 MacDonald Dissector, double ended, 190mm


13940-25 MacDonald Dissector, double ended, 250mm




13950-12 Watson Cheyne Dissector, fine with probe, 125mm

13950-18 Watson Cheyne Dissector, fine with probe, 180mm

Dressing Forceps



14000-12 Bryant Dressing Forcep, serrated jaws, box joint, 125mm




14020-12 French Bow Dressing Forcep, screw joint, 125mm



HEALTH MISSION
INSTRUMENTS


General Hernia Forceps

- 14300-14 Collingwood Stewart Hernia Forcep, ring bladed, screw joint, 140mm




Instrument Pins

- 14400-10 Mayo Instrument Pin, with ball end, 100mm
14400-15 Mayo Instrument Pin, with ball end, 150mm




Intestinal Forceps

- 14500-19 Denis Brown Intestinal Forcep, 190mm



Introducing Forceps

- 14600-19 Magill Introducing Forcep, screw joint, child, 190mm
14600-23 Magill Introducing Forcep, screw joint, adult, 230mm



General Knives



- 14750-20 Brain Knife, solid forged, 200mm blade
14750-22 Brain Knife, solid forged, 225mm blade
14750-25 Brain Knife, solid forged, 250mm blade
14750-30 Brain Knife, solid forged, 300mm blade



- 14760-20 Cartilage Knife, solid forged, 200mm, 70mm blade
14760-21 Cartilage Knife, solid forged, 215mm, 90mm blade



- 14770-10 Liston Amputation Knife, solid forged, 100mm blade
14770-12 Liston Amputation Knife, solid forged, 125mm blade
14770-15 Liston Amputation Knife, solid forged, 150mm blade
14770-18 Liston Amputation Knife, solid forged, 180mm blade




- 14780-14 Paget Knife, 140mm

Measures



- 15000-15 Stainless Steel Rule, 150mm
15000-30 Stainless Steel Rule, 300mm

Mouth Gags



- 15110-17 Ferguson Mouth Gag, Ackland jaw, child, 170mm
15110-20 Ferguson Mouth Gag, Ackland jaw, adult, 200mm




- 15120-17 Ferguson Mouth Gag, ordinary jaw, child, 170mm
15120-20 Ferguson Mouth Gag, ordinary jaw, adult, 200mm




General Needle Holders


- 15200-20 Bozemann Needleholder, double curve, serrated jaws, 200mm




- 15210-15 Crile Wood Needleholder, serrated jaws, 150mm
15210-18 Crile Wood Needleholder, serrated jaws, 180mm
15210-20 Crile Wood Needleholder, serrated jaws, 200mm
15210-23 Crile Wood Needleholder, serrated jaws, 230mm
15210-25 Crile Wood Needleholder, serrated jaws, 255mm
15210-30 Crile Wood Needleholder, serrated jaws, 305mm




- 15220-14 Crile-Murray Needleholder, serrated jaws, 145mm




- 15230-12 Derf Needleholder, serrated jaws, 120mm



General Needle Holders



15240-12 Halsey Needleholder, serrated jaws, 125mm




15250-20 Hegar Needleholder, serrated jaws, 200mm



15260-13 Kilner Needleholder, serrated jaws, 135mm

15260-16 Kilner Needleholder, serrated jaws, 165mm




15270-15 Lawrance Needleholder, lightweight, serrated jaws, 150mm




General Needle Holders


- 15280-15 Macphail Needleholder, screw joint with copper jaws, 150mm
15280-18 Macphail Needleholder, screw joint with copper jaws, 180mm
15280-20 Macphail Needleholder, screw joint with copper jaws, 200mm




- 15290-17 Mathieu Needleholder, serrated jaws, 170mm
15290-20 Mathieu Needleholder, serrated jaws, 200mm




- 15300-15 Mayo Needleholder, light pattern, serrated jaws, 150mm



- 15310-12 Mayo Needleholder, wide serrated jaws, 125mm
15310-14 Mayo Needleholder, wide serrated jaws, 140mm
15310-16 Mayo Needleholder, wide serrated jaws, 160mm
15310-18 Mayo Needleholder, wide serrated jaws, 185mm
15310-20 Mayo Needleholder, wide serrated jaws, 200mm
15310-23 Mayo Needleholder, wide serrated jaws, 230mm
15310-25 Mayo Needleholder, wide serrated jaws, 255mm
15310-30 Mayo Needleholder, wide serrated jaws, 305mm




General Needle Holders




15320-18 Nauntun Morgan Needleholder, compound action, serrated jaws, 185mm

15320-25 Nauntun Morgan Needleholder, compound action, serrated jaws, 250mm




15330-20 Thomson-Walker Needleholder, serrated jaws, 200mm



15340-18 Ward Needleholder, serrated, fenestrated jaws, 180mm

Needleholders - Tungsten Carbide



15400-15 Crile Wood Needleholder, tungsten carbide, 150mm

15400-18 Crile Wood Needleholder, tungsten carbide, 180mm

15400-20 Crile Wood Needleholder, tungsten carbide, 200mm

15400-23 Crile Wood Needleholder, tungsten carbide, 230mm

15400-25 Crile Wood Needleholder, tungsten carbide, 255mm

15400-30 Crile Wood Needleholder, tungsten carbide, 305mm



General Needle Holders-Tungsten Carbide


15410-18 DeBakey Needleholder, tungsten carbide, 180mm

15410-20 DeBakey Needleholder, tungsten carbide, 200mm


15410-23 DeBakey Needleholder, tungsten carbide, 230mm

15410-25 DeBakey Needleholder, tungsten carbide, 250mm

15410-30 DeBakey Needleholder, tungsten carbide, 305mm




15420-12 Derf Needleholder, tungsten carbide, 120mm




15430-20 Finochietto Needleholder, tungsten carbide, 200mm


15430-26 Finochietto Needleholder, tungsten carbide, 265mm



15440-12 Halsey Needleholder, tungsten carbide, 125mm




General Needle Holders-Tungsten Carbide




15460-13 Kilner Needleholder, tungsten carbide, 135mm

15460-16 Kilner Needleholder, tungsten carbide, 165mm



15470-15 Lawrance Needleholder, tungsten carbide, 150mm



15490-12 Mayo Needleholder, tungsten carbide, 125mm

15490-14 Mayo Needleholder, tungsten carbide, 140mm

15490-16 Mayo Needleholder, tungsten carbide, 160mm


15490-18 Mayo Needleholder, tungsten carbide, 185mm

15490-20 Mayo Needleholder, tungsten carbide, 200mm

15490-23 Mayo Needleholder, tungsten carbide, 230mm

15490-25 Mayo Needleholder, tungsten carbide, 255mm

15490-30 Mayo Needleholder, tungsten carbide, 305mm




15495-12 Neivert Needleholder, tungsten carbide, 125mm



General Needle Holders-Tungsten Carbide

- 15500-12 Ryder Needleholder, french eye, tungsten carbide, 125mm
- 15500-15 Ryder Needleholder, french eye, tungsten carbide, 150mm
- 15500-18 Ryder Needleholder, french eye, tungsten carbide, 180mm
- 15500-20 Ryder Needleholder, french eye, tungsten carbide, 200mm
- 15500-23 Ryder Needleholder, french eye, tungsten carbide, 230mm
- 15500-25 Ryder Needleholder, french eye, tungsten carbide, 255mm



Needles

- 15600-14 Aneurysm Needle, small curve, 145mm




- 15610-21 Bonney Reverdin Needle, handle only, standard model, without needles, 210mm



- 15620-20 Bonney Reverdin Needle, blunt, medium
- 15630-10 Bonney Reverdin Needle, sharp, large
- 15630-20 Bonney Reverdin Needle, sharp, medium
- 15630-30 Bonney Reverdin Needle, sharp, small



- 15640-19 Cooper Needle, 190mm




- 15650-27 Deschamps Needle, curved to left, 270mm




- 15660-27 Deschamps Needle, curved to right, 270mm



General Needles



- 15670-17 Dupuytren Needle, curved to left, 170mm
15680-17 Dupuytren Needle, curved to right, 170mm



- 15690-10 Jolis Aneurysm Needle, 200mm, No.1
15690-20 Jolis Aneurysm Needle, 200mm, No.2
15690-30 Jolis Aneurysm Needle, 200mm, No.3



- 15700-16 Symes Needle, 165mm

Percussors




- 15800-33 National Hospital Percussor, nylon handle, large rubber ring, 70mm diameter, 335mm

Plaster Instruments



- 15900-25 Bohler Bender, 255mm




- 15910-24 Bruns / Bohler Plaster Shears, saw edge, 240mm




General Plaster Instruments


15920-25 Guys Plaster Shears, flattened end, 250mm



15930-29 Guys Plaster Shears, probe end, 290mm




15940-28 Hennig Spreader, 280mm



15960-23 Stille-Lorenz Plaster Shears, 230mm

15960-26 Stille-Lorenz Plaster Shears, 260mm

15960-37 Stille-Lorenz Plaster Shears, 370mm



Probes

16030-35 Moynihan Gallstone Probe And Scoop, malleable, 350mm



HEALTH MISSION
INSTRUMENTS

General Probes



- 16040-12 Nickel Silver Probe With Eye, 125mm
- 16040-15 Nickel Silver Probe With Eye, 150mm
- 16040-18 Nickel Silver Probe With Eye, 180mm
- 16040-20 Nickel Silver Probe With Eye, 205mm
- 16040-23 Nickel Silver Probe With Eye, 230mm
- 16040-25 Nickel Silver Probe With Eye, 255mm
- 16040-30 Nickel Silver Probe With Eye, 305mm
- 16045-12 Silver Probe With Eye, 125mm
- 16045-15 Silver Probe With Eye, 150mm
- 16045-18 Silver Probe With Eye, 180mm
- 16045-20 Silver Probe With Eye, 205mm
- 16045-23 Silver Probe With Eye, 230mm
- 16045-25 Silver Probe With Eye, 255mm
- 16045-30 Silver Probe With Eye, 305mm




- 16050-15 Sinus Probe, silver, 150mm



- 16060-10 Stacks Probe, 105mm

Retractors




- 16100-24 Blake Amputation Flap Retractor, 240mm diameter



General Retractors


16105-70 Broomhead Psoas Retractor, 250mm, 22mm wide x 70mm deep blade

16105-90 Broomhead Psoas Retractor, 250mm, 28mm wide x 90mm deep blade



16110-10 Canny Ryall Army Pattern Retractor, large single blunt hook, small double blunt hook

16110-20 Canny Ryall Army Pattern Retractor, small single blunt hook, large double blunt hook



16120-15 Canny Ryall Navy Pattern Retractor, single hook 5mm wide, fenestrated blade 15mm wide

16120-20 Canny Ryall Navy Pattern Retractor, single hook 5mm wide, fenestrated blade 20mm wide

16120-25 Canny Ryall Navy Pattern Retractor, single hook 5mm wide, fenestrated blade 25mm wide

16120-30 Canny Ryall Navy Pattern Retractor, single hook 5mm wide, fenestrated blade 30mm wide



16130-18 Czerny Retractor, double ended, 180mm



16140-02 Deaver Abdominal Retractor, flat handle, No.00, 230mm, 25mm blade

16140-05 Deaver Abdominal Retractor, flat handle, No.0, 250mm, 25mm blade

16140-10 Deaver Abdominal Retractor, flat handle, No.1, 305mm, 25mm blade


16140-20 Deaver Abdominal Retractor, flat handle, No.2, 330mm, 25mm blade

16140-30 Deaver Abdominal Retractor, flat handle, No.3, 305mm, 38mm blade


16140-40 Deaver Abdominal Retractor, flat handle, No.4, 305mm, 50mm blade

16140-50 Deaver Abdominal Retractor, flat handle, No.5, 305mm, 75mm blade


16140-90 Deaver Abdominal Retractor, flat handle, 305mm, set of 5 (No.1 - No.5)




General Retractors



- 16150-02 Deaver Abdominal Retractor, hollow handle, No.00, 230mm, 25mm blade
- 16150-05 Deaver Abdominal Retractor, hollow handle, No.0, 250mm, 25mm blade
- 16150-10 Deaver Abdominal Retractor, hollow handle, No.1, 305mm, 25mm blade
- 16150-20 Deaver Abdominal Retractor, hollow handle, No.2, 330mm, 25mm blade
- 16150-30 Deaver Abdominal Retractor, hollow handle, No.3, 305mm, 38mm blade
- 16150-40 Deaver Abdominal Retractor, hollow handle, No.4, 305mm, 50mm blade
- 16150-50 Deaver Abdominal Retractor, hollow handle, No.5, 305mm, 75mm blade
- 16150-90 Deaver Abdominal Retractor, hollow handle, 305mm, set of 5 (No.1 - No.5)



- 16160-09 Doyen Retractor, on hollow handle, 90mm wide blade
- 16160-10 Doyen Retractor, on hollow handle, 100mm wide blade
- 16160-11 Doyen Retractor, on hollow handle, 115mm wide blade



- 16170-09 Doyen Retractor, on solid handle, 90mm wide blade
- 16170-10 Doyen Retractor, on solid handle, 100mm wide blade
- 16170-11 Doyen Retractor, on solid handle, 115mm wide blade




- 16180-10 Durham Retractor, 10mm blade
- 16180-16 Durham Retractor, 16mm blade
- 16180-22 Durham Retractor, 22mm blade




General Retractors


- 16190-24** Dyball And Lock Retractor, single ended, 240mm,
63mm x 112mm blade




- 16195-34** Dyson Retractor, 34mm wide



- 16205-23** Jolls Thyroid Retractor, 30mm x 14mm, 230mm
16206-23 Jolls Thyroid Retractor, 40mm x 20mm, 230mm




General Retractors




16210-38 Kelly Retractor, 270mm overall length, 190mm x 38mm blade

16210-57 Kelly Retractor, 270mm overall length, 190mm x 57mm blade




16211-38 Kelly Retractor, thackray pattern, 270mm overall length, 190mm x 38mm blade

16211-57 Kelly Retractor, thackray pattern, 270mm overall length, 190mm x 57mm blade



16220-11 Kocher Thyroid Retractor, shallow blade, 11mm wide



16220-14 Kocher Thyroid Retractor, deep blade, 14mm wide




16230-13 Lane Retractor, double ended, 230mm, 13mm x 19mm blade

16230-25 Lane Retractor, double ended, 230mm, 25mm x 32mm blade



General Retractors


- 16240-06 Langenbeck Retractor, single ended, 22mm x 6mm blade
- 16240-08 Langenbeck Retractor, single ended, 25mm x 8mm blade
- 16240-11 Langenbeck Retractor, single ended, 35mm x 11mm blade
- 16240-13 Langenbeck Retractor, single ended, 45mm x 13mm blade
- 16240-15 Langenbeck Retractor, single ended, 35mm x 15mm blade
- 16240-20 Langenbeck Retractor, single ended, 90mm x 20mm blade
- 16240-25 Langenbeck Retractor, single ended, 65mm x 25mm blade




- 16250-30 Morris Retractor, double ended, 215mm, 30mm & 40mm wide x 50mm deep blades
- 16250-40 Morris Retractor, double ended, 215mm, 40mm & 65mm wide x 50mm deep blades
- 16250-45 Morris Retractor, double ended, 215mm, 45mm & 50mm wide x 50mm deep blades
- 16255-30 Morris Retractor, double ended, 215mm, 30mm & 40mm wide x 70mm deep blades
- 16255-45 Morris Retractor, double ended, 215mm, 45mm & 50mm wide x 70mm deep blades
- 16255-60 Morris Retractor, double ended, 215mm, 60mm & 65mm wide x 70mm deep blades



- 16260-40 Morris Retractor, single ended, 250mm, 40mm wide x 50mm deep blade
- 16260-50 Morris Retractor, single ended, 250mm, 50mm wide x 50mm deep blade
- 16260-65 Morris Retractor, single ended, 250mm, 65mm wide x 50mm deep blade
- 16270-40 Morris Retractor, single ended, 250mm, 40mm wide x 70mm deep blade
- 16270-50 Morris Retractor, single ended, 250mm, 50mm wide x 70mm deep blade
- 16270-65 Morris Retractor, single ended, 250mm, 65mm wide x 70mm deep blade




- 16271-30 Morris Toothed Retractor, single ended, 215mm, 30mm wide x 80mm deep blade




General Retractors



16280-13 Moynihan Retractor, single ended, 215mm, 13mm x 22mm blade



16283-20 Ogilvie Retractor, 75mm expansion, 20mm wide x 20mm deep blades



16285-40 Oldfield Retractor, 4 prong, 270mm, broad blade, 50mm wide


16285-50 Oldfield Retractor, 5 prong, 270mm, narrow blade, 40mm wide



16295-13 Ollier Retractor, single ended, 2 prong, 240mm, 13mm wide x 20mm deep blade

16295-25 Ollier Retractor, single ended, 2 prong, 240mm, 25mm wide x 45mm deep blade

16295-37 Ollier Retractor, single ended, 2 prong, 240mm, 37mm wide x 60mm deep blade



16320-21 Parker-Langenbeck (U S Army) Retractor, set of 2, 210mm



16330-31 Paton Cartilage Retractor, 31mm x 10mm


16330-44 Paton Cartilage Retractor, 44mm x 10mm

16330-51 Paton Cartilage Retractor, 51mm x 14mm



General Retractors

16332-65 Proust Retractor, 85mm x 45mm




16334-10 Pyrah Retractor, 75mm x 50mm

16334-20 Pyrah Retractor, 100mm x 50mm

16334-30 Pyrah Retractor, 100mm x 65mm

16334-40 Pyrah Retractor, 125mm x 70mm

16334-50 Pyrah Retractor, 150mm x 80mm



16336-03 Shucksmith Retractor, 30mm x 13mm

16336-05 Shucksmith Retractor, 50mm x 13mm

16336-06 Shucksmith Retractor, 65mm x 30mm

16336-08 Shucksmith Retractor, 80mm x 30mm

16336-09 Shucksmith Retractor, 95mm x 45mm

16336-11 Shucksmith Retractor, 115mm x 50mm


16336-12 Shucksmith Retractor, 125mm x 55mm

16336-14 Shucksmith Retractor, 140mm x 60mm


16336-15 Shucksmith Retractor, 155mm x 70mm

16336-17 Shucksmith Retractor, 170mm x 80mm


16336-90 Shucksmith Retractor, set of 10 complete with storage clamp




General Retractors



- 16340-12 St. Marks Hospital Retractor, with lip, 125mm x 60mm blade
- 16340-17 St. Marks Hospital Retractor, with lip, 175mm x 60mm blade
- 16342-12 St. Marks Hospital Retractor, with lip, 125mm x 60mm blade, with twin fibrelight
- 16342-17 St. Marks Hospital Retractor, with lip, 175mm x 60mm blade, with twin fibrelight



- 16340-20 St. Marks Hospital Retractor, with lip, 205mm x 60mm blade, cranked
- 16342-20 St. Marks Hospital Retractor, with lip, 205mm x 60mm blade, cranked, with twin fibrelight
- 16345-20 St. Marks Hospital Retractor, twin fibrelight stem only



- 16350-12 St. Marks Hospital Retractor, without lip, 125mm x 60mm blade
- 16350-17 St. Marks Hospital Retractor, without lip, 175mm x 60mm blade
- 16352-12 St. Marks Hospital Retractor, without lip, 125mm x 60mm blade, with twin fibrelight
- 16352-17 St. Marks Hospital Retractor, without lip, 175mm x 60mm blade, with twin fibrelight



- 16358-13 Stack Retractor, 13mm deep blades, 140mm
- 16358-25 Stack Retractor, 25mm deep blades, 140mm



General Retractors


- 16360-18 Tracheal Hook Retractor, double, blunt, 180mm
16370-18 Tracheal Hook Retractor, double, sharp, 180mm




- 16380-18 Tracheal Hook Retractor, single, blunt, 150mm
16390-18 Tracheal Hook Retractor, single, sharp, 150mm



- 16400-10 Volkmann Retractor, blunt, 220mm, 1 prong
16400-20 Volkmann Retractor, blunt, 220mm, 2 prong
16400-30 Volkmann Retractor, blunt, 220mm, 3 prong
16400-40 Volkmann Retractor, blunt, 220mm, 4 prong
16400-50 Volkmann Retractor, blunt, 220mm, 5 prong
16400-60 Volkmann Retractor, blunt, 220mm, 6 prong
16410-10 Volkmann Retractor, semi-sharp, 220mm, 1 prong
16410-20 Volkmann Retractor, semi-sharp, 220mm, 2 prong
16410-30 Volkmann Retractor, semi-sharp, 220mm, 3 prong
16410-40 Volkmann Retractor, semi-sharp, 220mm, 4 prong
16410-50 Volkmann Retractor, semi-sharp, 220mm, 5 prong
16410-60 Volkmann Retractor, semi-sharp, 220mm, 6 prong
16420-10 Volkmann Retractor, sharp, 220mm, 1 prong
16420-20 Volkmann Retractor, sharp, 220mm, 2 prong
16420-30 Volkmann Retractor, sharp, 220mm, 3 prong
16420-40 Volkmann Retractor, sharp, 220mm, 4 prong
16420-50 Volkmann Retractor, sharp, 220mm, 5 prong
16420-60 Volkmann Retractor, sharp, 220mm, 6 prong






16430-06	Walton Retractor, copper, malleable, 180mm x 6mm
16430-08	Walton Retractor, copper, malleable, 180mm x 8mm
16430-11	Walton Retractor, copper, malleable, 180mm x 11mm
16430-12	Walton Retractor, copper, malleable, 180mm x 12mm
16430-13	Walton Retractor, copper, malleable, 180mm x 13mm
16430-16	Walton Retractor, copper, malleable, 180mm x 16mm
16430-19	Walton Retractor, copper, malleable, 180mm x 19mm
16430-25	Walton Retractor, copper, malleable, 180mm x 25mm
16430-32	Walton Retractor, copper, malleable, 180mm x 32mm
16430-38	Walton Retractor, copper, malleable, 180mm x 38mm
16430-50	Walton Retractor, copper, malleable, 180mm x 50mm
16430-63	Walton Retractor, copper, malleable, 180mm x 63mm
16435-08	Walton Retractor, copper, malleable, 200mm x 8mm
16435-13	Walton Retractor, copper, malleable, 200mm x 13mm
16435-16	Walton Retractor, copper, malleable, 200mm x 16mm
16435-19	Walton Retractor, copper, malleable, 200mm x 19mm
16435-25	Walton Retractor, copper, malleable, 200mm x 25mm
16440-13	Walton Retractor, copper, malleable, 255mm x 13mm
16440-19	Walton Retractor, copper, malleable, 255mm x 19mm
16440-25	Walton Retractor, copper, malleable, 255mm x 25mm
16440-32	Walton Retractor, copper, malleable, 255mm x 32mm
16440-38	Walton Retractor, copper, malleable, 255mm x 38mm
16440-50	Walton Retractor, copper, malleable, 255mm x 50mm
16440-63	Walton Retractor, copper, malleable, 255mm x 63mm
16450-13	Walton Retractor, copper, malleable, 300mm x 13mm
16450-19	Walton Retractor, copper, malleable, 300mm x 19mm
16450-25	Walton Retractor, copper, malleable, 300mm x 25mm
16450-32	Walton Retractor, copper, malleable, 300mm x 32mm
16450-38	Walton Retractor, copper, malleable, 300mm x 38mm
16450-45	Walton Retractor, copper, malleable, 300mm x 45mm
16450-50	Walton Retractor, copper, malleable, 300mm x 50mm
16450-63	Walton Retractor, copper, malleable, 300mm x 63mm




General Retractors

- 16460-50** Walton Retractor, tapered, copper, malleable.
300mm x 50mm tapering to 25mm





Retractors - Self Retaining

- 16500-09** Balfour Self Retaining Retractor, fenestrated end blades, baby, 90mm
16500-10 Balfour Self Retaining Retractor, fenestrated end blades, child, 100mm
16500-18 Balfour Self Retaining Retractor, fenestrated end blades, adult, 180mm




- 16505-19** Balfour Self Retaining Retractor, with trigger ratchet, 2 pairs detachable end blades 45mm wide x 65mm & 85mm deep, centre blade 65mm x 65mm, 195mm





16510-60 Berkeley Abdominal Self Retaining Retractor, parallel action, 35mm, 60mm blades



16520-10 Cleland Self Retaining Retractor, deep angled blade only, 115mm x 65mm

16520-20 Cleland Self Retaining Retractor, frame only, 250mm ID

16520-30 Cleland Self Retaining Retractor, shallow curved blade only, 50mm x 40mm


16520-90 Cleland Self Retaining Retractor, complete set, 1 x frame, 4 x shallow curved blades, 4 x deep angled blades




General Retractors-Self Retaining

16530-90 Collin Self Retaining Retractor, complete with centre blade, 38mm, 60mm blades

16540-90 Collin Self Retaining Retractor, without centre blade, 38mm, 60mm blades




16545-85 DeBakey-Balfour Self Retaining Retractor, with trigger ratchet, blades 45mm wide x 85mm deep, centre blade 85mm wide x 70mm deep




16550-50 Denis Browne Self Retaining Retractor, oval frame only, 180mm x 140mm

16550-90 Denis Browne Self Retaining Retractor, set, 1 x frame, 4 x blades A B C, 2 x blade D, 1 x blades E F G H J




HEALTH MISSION
INSTRUMENTS


General Retractors-Self Retaining




16550-05 Denis Browne Self Retaining Retractor, blade 'A'
only, curved claw blade, 13mm x 13mm




16550-10 Denis Browne Self Retaining Retractor, blade 'B'
only, angled claw blade, 60mm x 20mm



16550-15 Denis Browne Self Retaining Retractor, blade 'C'
only, curved blade, 40mm x 35mm



16550-20 Denis Browne Self Retaining Retractor, blade 'D'
only, curved hook blade, sharp



16550-25 Denis Browne Self Retaining Retractor, blade 'E'
only, lung retractor end to fit shaft 'H'



General Retractors-Self Retaining

- 16550-30 Denis Browne Self Retaining Retractor, blade 'F' only, curved blade end to fit shaft 'H'




- 16550-35 Denis Browne Self Retaining Retractor, blade 'G' only, regressive curved blade end to fit shaft 'H'




- 16550-40 Denis Browne Self Retaining Retractor, shaft 'H' only, with threaded retaining action




- 16550-45 Denis Browne Self Retaining Retractor, spanner 'J' only, to fit shaft 'H'




General Retractors-Self Retaining



- 16560-10 Goligher Self Retaining Retractor, wing nut only
- 16560-20 Goligher Self Retaining Retractor, knurled end nut only
- 16560-50 Goligher Self Retaining Retractor, deep lateral blade only, 50mm wide x 90mm deep
- 16560-55 Goligher Self Retaining Retractor, shallow lateral blade only, 50mm
- 16560-60 Goligher Self Retaining Retractor, deep centre blade only, 90mm wide x 90mm deep
- 16560-65 Goligher Self Retaining Retractor, shallow centre blade only, 90mm wide x 55mm deep
- 16560-70 Goligher Self Retaining Retractor, wingnut and mount for centre blade only
- 16560-80 Goligher Self Retaining Retractor, frame only
- 16560-90 Goligher Self Retaining Retractor, complete with 2 centre blades, 2 pairs lateral blades



- 16570-90 Gossett Panting Self Retaining Retractor, large, complete with centre blade
- 16580-90 Gossett Panting Self Retaining Retractor, small, complete with centre blade




- 16590-90 Gossett Self Retaining Retractor, large
- 16600-90 Gossett Self Retaining Retractor, small




General Retractors-Self Retaining


- 16605-18** Irwin Self Retaining Retractor, cam rack, 4/4 teeth, 180mm




- 16610-15** Jolls Thyroid Self Retaining Retractor, 150mm




- 16620-90** Lace Self Retaining Retractor, with 1 pair side blades




- 16625-90** McKeown Self Retaining Retractor, 60mm deep x 45mm wide blades






16630-21 Norfolk And Norwich Self Retaining Retractor, cam rack, 4/5 teeth, 215mm, 50mm blades




16640-20 Ricard Begouin Self Retaining Retractor, 200mm spread, 35mm deep x 65mm wide and 65mm deep x 65mm wide blades

16645-28 Ricard Begouin Griffiths Self Retaining Retractor, 280mm spread, 60mm deep x 80mm wide and 100mm deep x 80mm wide blades



16650-05 Swift Joly Self Retaining Retractor, blade 'A' only, 50mm deep x 25mm wide




16650-10 Swift Joly Self Retaining Retractor, blade 'B' only, 70mm deep x 25mm wide




General Retractors-Self Retaining

- 16650-15** Swift Joly Self Retaining Retractor, blade 'C' only,
90mm deep x 25mm wide




- 16650-20** Swift Joly Self Retaining Retractor, blade 'D' only,
65mm deep x 65mm wide




- 16650-50** Swift Joly Self Retaining Retractor, oval frame
only, 180mm x 140mm


- 16650-90** Swift Joly Self Retaining Retractor, complete, 1 x
frame, 1 x blades A B, 2 x blades C D



General Retractors-Self Retaining



- 16655-90 Stevens Abdominal Self Retaining Retractor, lightweight, 60mm deep x 45mm wide blades, 190g (7oz)



- 16658-05 St. Thomas Hospital Self Retaining Retractor, blade 'A' only, right angled, 230mm long with 9 pegs, 145mm deep x 50mm wide
- 16658-10 St. Thomas Hospital Self Retaining Retractor, blade 'B' only, right angled, 215mm long with 9 pegs, 125mm deep x 38mm wide
- 16658-15 St. Thomas Hospital Self Retaining Retractor, blade 'C' only, right angled, 180mm long with 6 pegs, 120mm deep x 38mm wide
- 16658-20 St. Thomas Hospital Self Retaining Retractor, blade 'D' only, right angled, 180mm long with 6 pegs, 95mm deep x 38mm wide




- 16658-25 St. Thomas Hospital Self Retaining Retractor, blade 'E' only, curved, 215mm long with 9 pegs, 45mm deep x 38mm wide
- 16658-30 St. Thomas Hospital Self Retaining Retractor, blade 'F' only, curved, 180mm long with 7 pegs, 45mm deep x 38mm wide
- 16658-35 St. Thomas Hospital Self Retaining Retractor, blade 'G' only, curved, 125mm long with 4 pegs, 45mm deep x 38mm wide
- 16658-50 St. Thomas Hospital Self Retaining Retractor, blade 'K' only, curved, 215mm long with 9 pegs, 65mm deep x 38mm wide
- 16658-52 St. Thomas Hospital Self Retaining Retractor, blade 'L' only, curved, 180mm long with 7 pegs, 65mm deep x 38mm wide
- 16658-54 St. Thomas Hospital Self Retaining Retractor, blade 'M' only, curved, 125mm long with 4 pegs, 65mm deep x 38mm wide
- 16658-56 St. Thomas Hospital Self Retaining Retractor, blade 'N' only, curved, 215mm long with 9 pegs, 85mm deep x 38mm wide



General Retractors-Self Retaining

16658-40 St.Thomas Hospital Self Retaining Retractor, blade 'H' only, deaver type, guttered, 230mm long with 6 pegs, 115mm deep x 70mm tapering to 64mm


16658-45 St.Thomas Hospital Self Retaining Retractor, blade 'J' only, deaver type, guttered, 255mm long with 6 pegs, 150mm deep x 70mm tapering to 64mm




16658-60 St.Thomas Hospital Self Retaining Retractor, frame 'A' only, 255mm x 345mm

16658-70 St.Thomas Hospital Self Retaining Retractor, frame 'B' only, 330mm x 445mm


16658-90 St.Thomas Hospital Self Retaining Retractor, complete, 1 x frame A, 1 x blades A B C D E F G H J




16660-20 Travers Self Retaining Retractor, 200mm




16670-05 Turner Warwick Universal Self Retaining Retractor, blade 'A' only, copper blade, 19mm x 76mm




General Retractors-Self Retaining



16670-10 Turner Warwick Universal Self Retaining Retractor, blade 'B' only, copper blade, with clip, 38mm x 102mm



16670-15 Turner Warwick Universal Self Retaining Retractor, blade 'C' only, copper blade, with clip, 57mm x 89mm



16670-20 Turner Warwick Universal Self Retaining Retractor, blade 'D' only, copper blade, with clip, 57mm x 114mm




General Retractors-Self Retaining


- 16670-25** Turner Warwick Universal Self Retaining Retractor, blade 'E' only, curved, 51mm dia. x 38mm




- 16670-30** Turner Warwick Universal Self Retaining Retractor, blade 'F' only, curved, with clip, 76mm dia. x 60mm



- 16670-35** Turner Warwick Universal Self Retaining Retractor, blade 'G' only, flat, with clip, 89mm x 60mm




- 16670-40** Turner Warwick Universal Self Retaining Retractor, blade 'H' only, flat, with clip, 114mm x 60mm




General Retractors-Self Retaining




16670-45 Turner Warwick Universal Self Retaining Retractor, blade 'J' only, single hook, 13mm deep



16670-50 Turner Warwick Universal Self Retaining Retractor, blade 'K' only, flat, with clip, 55mm x 90mm



16670-52 Turner Warwick Universal Self Retaining Retractor, blade 'L' only, flat, with clip, 55mm x 120mm




16670-54 Turner Warwick Universal Self Retaining Retractor, blade 'M' only, copper blade, 31mm x 55mm




General Retractors-Self Retaining


- 16670-56** Turner Warwick Universal Self Retaining Retractor, blade 'N' only, deep curved, 60mm x 55mm




- 16670-58** Turner Warwick Universal Self Retaining Retractor, blade 'P' only, shallow curved, 28mm x 30mm




- 16670-60** Turner Warwick Universal Self Retaining Retractor, circular frame 'A' only, 178mm ID
16670-65 Turner Warwick Universal Self Retaining Retractor, circular frame 'B' only, 229mm ID
16670-70 Turner Warwick Universal Self Retaining Retractor, circular frame 'C' only, 305mm ID
16670-75 Turner Warwick Universal Self Retaining Retractor, oval frame 'D' only, 457mm x 330mm ID
16670-80 Turner Warwick Universal Self Retaining Retractor, oval frame 'E' only, 370mm x 300mm ID
16670-90 Turner Warwick Universal Self Retaining Retractor, complete, 1 x frames A B C D, 2 x blades A B C D E F G H J
16670-95 Turner Warwick Universal Self Retaining Retractor, complete, 1 x frames B C, 4 x blade E, 2 x blade N, 1 x blades C K L



- 16680-14** Weislander Self Retaining Retractor, 3/4 teeth, blunt, 140mm
16690-14 Weislander Self Retaining Retractor, 3/4 teeth, sharp, 140mm




General Retractors-Self Retaining




16700-14 West Self Retaining Retractor, sharp, 140mm

Saws



16900-20 Amputation Saw, 200mm

16900-22 Amputation Saw, 225mm



16920-10 Engels Plaster Saw

Scalpel Handles




17020-30 Scalpel Handle, No.3



17020-35 Scalpel Handle, No.3L



17020-40 Scalpel Handle, No.4



17020-45 Scalpel Handle, No.4L



HEALTH MISSION
INSTRUMENTS

General Scalpel Handles

17020-70 Scalpel Handle, No.7



Scissors


17100-15 Aufricht Dissecting Scissor, curved, light pattern, 150mm

17120-15 Aufricht Dissecting Scissor, straight, light pattern, 150mm




17110-14 Aufricht Dissecting Scissor, curved, triangular blade, 140mm

17130-14 Aufricht Dissecting Scissor, straight, triangular blade, 140mm




17140-16 Carless Scissor, straight, lightweight, blunt point, 165mm

17140-20 Carless Scissor, straight, lightweight, blunt point, 200mm



17160-18 Doyen Scissor, straight, blunt point, 180mm



17170-10 Dressing Scissor, curved, blunt / blunt, 100mm


17170-12 Dressing Scissor, curved, blunt / blunt, 125mm

17170-14 Dressing Scissor, curved, blunt / blunt, 140mm


17170-15 Dressing Scissor, curved, blunt / blunt, 150mm

17170-18 Dressing Scissor, curved, blunt / blunt, 180mm


17170-20 Dressing Scissor, curved, blunt / blunt, 200mm




General Scissors




- 17180-10 Dressing Scissor, curved, blunt / sharp, 100mm
- 17180-12 Dressing Scissor, curved, blunt / sharp, 125mm
- 17180-14 Dressing Scissor, curved, blunt / sharp, 140mm
- 17180-15 Dressing Scissor, curved, blunt / sharp, 150mm
- 17180-18 Dressing Scissor, curved, blunt / sharp, 180mm
- 17180-20 Dressing Scissor, curved, blunt / sharp, 200mm




- 17190-10 Dressing Scissor, curved, sharp / sharp, 100mm
- 17190-12 Dressing Scissor, curved, sharp / sharp, 125mm
- 17190-14 Dressing Scissor, curved, sharp / sharp, 140mm
- 17190-15 Dressing Scissor, curved, sharp / sharp, 150mm
- 17190-18 Dressing Scissor, curved, sharp / sharp, 180mm
- 17190-20 Dressing Scissor, curved, sharp / sharp, 200mm



- 17200-10 Dressing Scissor, straight, blunt / blunt, 100mm
- 17200-12 Dressing Scissor, straight, blunt / blunt, 125mm
- 17200-14 Dressing Scissor, straight, blunt / blunt, 140mm
- 17200-15 Dressing Scissor, straight, blunt / blunt, 150mm
- 17200-18 Dressing Scissor, straight, blunt / blunt, 180mm
- 17200-20 Dressing Scissor, straight, blunt / blunt, 200mm



- 17210-10 Dressing Scissor, straight, blunt / sharp, 100mm
- 17210-12 Dressing Scissor, straight, blunt / sharp, 125mm
- 17210-14 Dressing Scissor, straight, blunt / sharp, 140mm
- 17210-15 Dressing Scissor, straight, blunt / sharp, 150mm
- 17210-18 Dressing Scissor, straight, blunt / sharp, 180mm
- 17210-20 Dressing Scissor, straight, blunt / sharp, 200mm




- 17220-12 Dressing Scissor, straight, probe / sharp, 125mm
- 17220-14 Dressing Scissor, straight, probe / sharp, 145mm




General Scissors


- 17230-10 Dressing Scissor, straight, sharp / sharp, 100mm
- 17230-12 Dressing Scissor, straight, sharp / sharp, 125mm
- 17230-14 Dressing Scissor, straight, sharp / sharp, 140mm
- 17230-15 Dressing Scissor, straight, sharp / sharp, 150mm
- 17230-18 Dressing Scissor, straight, sharp / sharp, 180mm
- 17230-20 Dressing Scissor, straight, sharp / sharp, 200mm




-
- 17235-08 Elliot Iris Scissor, sharp points, 85mm



-
- 17240-14 Fergusson Scissor, 145mm
 - 17240-15 Fergusson Scissor, 150mm
 - 17240-18 Fergusson Scissor, 180mm




-
- 17245-10 Glasgow Scissor, straight, fine points, 100mm




-
- 17260-15 Heaths Stitch Scissor, 150mm




General Scissors




- 17270-14 Listers Bandage Scissor, 140mm
 17270-18 Listers Bandage Scissor, 180mm
 17270-20 Listers Bandage Scissor, 200mm
 17270-23 Listers Bandage Scissor, 235mm




- 17280-18 Listers Bandage Scissor, with double finger bow, 185mm
 17280-21 Listers Bandage Scissor, with double finger bow, 210mm



- 17290-16 Listers Eccentric Scissor, 160mm
 17290-19 Listers Eccentric Scissor, 190mm



- 17300-12 Littauer Scissor, miniature, straight with hook end, 120mm




- 17310-12 Littauer Scissor, straight with hook end, 125mm
 17310-14 Littauer Scissor, straight with hook end, 140mm
 17310-18 Littauer Scissor, straight with hook end, 180mm
 17380-14 Metzenbaum Scissor, fine, curved, blunt / blunt, 145mm
 17380-18 Metzenbaum Scissor, fine, curved, blunt / blunt, 180mm
 17380-20 Metzenbaum Scissor, fine, curved, blunt / blunt, 200mm
 17380-23 Metzenbaum Scissor, fine, curved, blunt / blunt, 230mm




General Scissors


- 17330-14 Mayo Scissor, curved, bevelled blades, 140mm
 17330-16 Mayo Scissor, curved, bevelled blades, 165mm
 17330-19 Mayo Scissor, curved, bevelled blades, 190mm
 17330-21 Mayo Scissor, curved, bevelled blades, 215mm




- 17340-14 Mayo Scissor, curved, round blades, 140mm
 17340-17 Mayo Scissor, curved, round blades, 170mm




- 17350-14 Mayo Scissor, straight, bevelled blades, 140mm
 17350-16 Mayo Scissor, straight, bevelled blades, 165mm
 17350-19 Mayo Scissor, straight, bevelled blades, 190mm
 17350-21 Mayo Scissor, straight, bevelled blades, 215mm




- 17360-14 Mayo Scissor, straight, round blades, 140mm
 17360-17 Mayo Scissor, straight, round blades, 170mm




- 17370-14 Metzenbaum Scissor, curved, blunt / blunt, 145mm
 17370-18 Metzenbaum Scissor, curved, blunt / blunt, 180mm
 17370-20 Metzenbaum Scissor, curved, blunt / blunt, 200mm
 17370-23 Metzenbaum Scissor, curved, blunt / blunt, 230mm
 17370-25 Metzenbaum Scissor, curved, blunt / blunt, 250mm
 17370-28 Metzenbaum Scissor, curved, blunt / blunt, 280mm
 17370-30 Metzenbaum Scissor, curved, blunt / blunt, 300mm



- 17390-14 Metzenbaum Scissor, fine, curved, sharp / sharp, 145mm



General Scissors




17400-14 Metzenbaum Scissor, fine, straight, blunt / blunt, 145mm

17400-18 Metzenbaum Scissor, fine, straight, blunt / blunt, 180mm

17400-20 Metzenbaum Scissor, fine, straight, blunt / blunt, 200mm

17400-23 Metzenbaum Scissor, fine, straight, blunt / blunt, 230mm



17410-14 Metzenbaum Scissor, fine, straight, sharp / sharp, 145mm



17420-14 Metzenbaum Scissor, straight, blunt / blunt, 145mm

17420-18 Metzenbaum Scissor, straight, blunt / blunt, 180mm

17420-20 Metzenbaum Scissor, straight, blunt / blunt, 200mm

17420-23 Metzenbaum Scissor, straight, blunt / blunt, 230mm

17420-25 Metzenbaum Scissor, straight, blunt / blunt, 250mm

17420-28 Metzenbaum Scissor, straight, blunt / blunt, 280mm

17420-30 Metzenbaum Scissor, straight, blunt / blunt, 300mm



17430-23 Nelson Lobectomy Scissor, curved, 230mm




17440-23 Nelson Lobectomy Scissor, straight, 230mm




General Scissors


- 17450-14 Post Mortem Scissor, ball end, open shank,
140mm




- 17460-10 Post Mortem Scissor, probe / sharp, open shank,
100mm




- 17490-10 Spencer Scissor, hook end, 105mm




- 17493-11 Stevens Scissor, curved, blunt, 115mm




- 17494-11 Stevens Scissor, curved, sharp, 115mm







17495-11 Stevens Scissor, straight, blunt, 115mm



17496-11 Stevens Scissor, straight, sharp, 115mm



17500-17 Universal Bandage Scissor, saw edge,
autoclavable handles, 175mm




17510-12 Universal Wire Cutting Scissor, saw edge, 120mm




General Scissors-Tungsten Carbide


- 17700-14 Dressing Scissor, tungsten carbide, curved, blunt / blunt, 145mm
 17710-14 Dressing Scissor, tungsten carbide, curved, blunt / sharp, 145mm
 17720-14 Dressing Scissor, tungsten carbide, curved, sharp / sharp, 145mm




- 17730-14 Dressing Scissor, tungsten carbide, straight, blunt / blunt, 145mm
 17740-14 Dressing Scissor, tungsten carbide, straight, blunt / sharp, 145mm
 17750-14 Dressing Scissor, tungsten carbide, straight, sharp / sharp, 145mm




- 17755-11 Iris Scissor, tungsten carbide, curved, sharp, 115mm
 17756-11 Iris Scissor, tungsten carbide, straight, sharp, 115mm



- 17760-14 Mayo Scissor, tungsten carbide, curved, 145mm
 17760-17 Mayo Scissor, tungsten carbide, curved, 170mm
 17760-19 Mayo Scissor, tungsten carbide, curved, 190mm
 17770-14 Mayo Scissor, tungsten carbide, straight, 145mm
 17770-17 Mayo Scissor, tungsten carbide, straight, 170mm
 17770-19 Mayo Scissor, tungsten carbide, straight, 190mm



General Scissors-Tungsten Carbide




17780-15 Mayo-Stille Scissor , tungsten carbide, curved, 150mm

17780-17 Mayo-Stille Scissor , tungsten carbide, curved, 170mm

17790-15 Mayo-Stille Scissor , tungsten carbide, straight, 150mm

17790-17 Mayo-Stille Scissor , tungsten carbide, straight, 170mm




17800-14 Metzenbaum Scissor, tungsten carbide, curved, blunt / blunt, 145mm

17800-17 Metzenbaum Scissor, tungsten carbide, curved, blunt / blunt, 180mm

17800-20 Metzenbaum Scissor, tungsten carbide, curved, blunt / blunt, 200mm

17800-23 Metzenbaum Scissor, tungsten carbide, curved, blunt / blunt, 230mm




17810-14 Metzenbaum Scissor, tungsten carbide, straight, blunt / blunt, 145mm


17810-17 Metzenbaum Scissor, tungsten carbide, straight, blunt / blunt, 180mm

17810-20 Metzenbaum Scissor, tungsten carbide, straight, blunt / blunt, 200mm

17810-23 Metzenbaum Scissor, tungsten carbide, straight, blunt / blunt, 230mm



17815-23 Nelson Lobectomy Scissor, tungsten carbide, curved, 230mm




17816-23 Nelson Lobectomy Scissor, tungsten carbide, straight, 230mm



General Scissors-Tungsten Carbide


- 17820-12 Universal Wire Cutting Scissor, tungsten carbide, 120mm



Scissors - Ultra Cutting


- 17600-14 Dressing Scissor, ultra cutting, curved, blunt / sharp, 145mm

- 17610-14 Dressing Scissor, ultra cutting, straight, blunt / sharp, 145mm




- 17620-11 Iris Scissor, ultra cutting, curved, 115mm

- 17630-11 Iris Scissor, ultra cutting, straight, 115mm




- 17635-11 Kilner Scissor, ultra cutting, curved, blunt, 110mm


- 17636-11 Kilner Scissor, ultra cutting, straight, sharp, 110mm




General Scissors-Ultra Cutting




- 17640-14 Mayo Scissor, ultra cutting, curved, 145mm
- 17640-17 Mayo Scissor, ultra cutting, curved, 170mm
- 17640-20 Mayo Scissor, ultra cutting, curved, 200mm
- 17640-21 Mayo Scissor, ultra cutting, curved, 215mm
- 17650-14 Mayo Scissor, ultra cutting, straight, 145mm
- 17650-17 Mayo Scissor, ultra cutting, straight, 170mm
- 17650-20 Mayo Scissor, ultra cutting, straight, 200mm
- 17650-21 Mayo Scissor, ultra cutting, straight, 215mm



- 17655-15 Mayo-Stille Scissor, ultra cutting, curved, 150mm
- 17655-17 Mayo-Stille Scissor, ultra cutting, curved, 170mm
- 17655-19 Mayo-Stille Scissor, ultra cutting, curved, 190mm
- 17656-15 Mayo-Stille Scissor, ultra cutting, straight, 150mm
- 17656-17 Mayo-Stille Scissor, ultra cutting, straight, 170mm
- 17656-19 Mayo-Stille Scissor, ultra cutting, straight, 190mm



- 17660-14 Metzenbaum Scissor, ultra cutting, curved, 145mm
- 17660-18 Metzenbaum Scissor, ultra cutting, curved, 180mm
- 17660-20 Metzenbaum Scissor, ultra cutting, curved, 200mm
- 17660-23 Metzenbaum Scissor, ultra cutting, curved, 230mm
- 17660-26 Metzenbaum Scissor, ultra cutting, curved, 260mm
- 17670-14 Metzenbaum Scissor, ultra cutting, straight, 145mm




- 17680-11 Stevens Scissor, ultra cutting, curved, 110mm
- 17681-11 Stevens Scissor, ultra cutting, straight, 110mm




General Scissors-Ultra Cutting

- 17690-11 Strabismus Scissor, ultra cutting, curved, 110mm
17691-11 Strabismus Scissor, ultra cutting, straight, 110mm




Sinus Forceps

- 18000-12 Lister Sinus Forcep, serrated tip, 125mm
18000-15 Lister Sinus Forcep, serrated tip, 150mm
18000-18 Lister Sinus Forcep, serrated tip, 180mm
18000-20 Lister Sinus Forcep, serrated tip, 200mm



Splinter Forceps

- 18200-12 Arthur Splinter Forcep, pointed, serrated jaws, 125mm




- 18210-10 Hunter Splinter Forcep, straight, 105mm



- 18220-07 Martin Splinter Forcep, 75mm
18220-09 Martin Splinter Forcep, 90mm
18220-11 Martin Splinter Forcep, 115mm




General Sponge Holding Forceps




18300-18 Foerster Sponge Holding Forcep, curved, serrated jaws, 180mm

18300-24 Foerster Sponge Holding Forcep, curved, serrated jaws, 240mm



18310-18 Foerster Sponge Holding Forcep, straight, serrated jaws, 180mm


18310-24 Foerster Sponge Holding Forcep, straight, serrated jaws, 240mm



18320-18 Rampley Sponge Holding Forcep, straight, serrated jaws, 180mm

18320-24 Rampley Sponge Holding Forcep, straight, serrated jaws, 240mm

Steriliser Forceps



18410-26 Cheattle Steriliser Forcep, light, box joint, 265mm

Suction Tubes



18600-10 American Pattern Suction Tube, s/s, with stillette, 165mm



18601-10 American Pattern Suction Tube, s/s, insulated, with stillette, 165mm




General Suction Tubes

18610-17 Lawrie Extension, s/s, 175mm



18620-10 Millin Suction Tube



18650-24 Poole Suction Tube, s/s, curved mount, 240mm, 10mm dia.




18660-24 Poole Suction Tube, s/s, straight mount, 240mm, 10mm dia.




18664-15 Poole Suction Tube, s/s, paediatric, straight mount, 155mm, 6mm dia.

18665-31 Simpson Smith Suction Tube, thumb button control, detachable inner suction tube and perforated sheath




18665-50 Simpson Smith Suction Tube, spare spring only

18670-25 Tinckler (Simpson Smith) Suction Tube, with on / off valve, 150mm, 255mm



General Suction Tubes




- 18680-10 Universal Suction Tube, handle only
- 18680-20 Universal Suction Tube, abdominal tube, fenestrated tip
- 18680-30 Universal Suction Tube, bent tube, oval open tip
- 18680-40 Universal Suction Tube, bent tube, round open tip
- 18680-50 Universal Suction Tube, malleable tube, perforated tip
- 18680-60 Universal Suction Tube, yankauer tube, perforated tip
- 18680-90 Universal Suction Tube, complete set, handle and 5 interchangeable ends




- 18700-30 Yankauer Suction Tube, s/s, 305mm
- 18700-50 Yankauer Suction Tube, s/s, spare tip only

Suture Forceps



- 18800-15 Guy's Suture Forcep, angled flat blades, screw joint, 150mm

Tissue Forceps




- 18900-30 Allis (American Pattern) Tissue Forcep, 160mm, 3/4 teeth




HEALTH MISSION
INSTRUMENTS

General Tissue Forceps


- 18910-20 Allis Tissue Forcep, 150mm, 2/3 teeth
- 18910-30 Allis Tissue Forcep, 150mm, 3/4 teeth
- 18910-40 Allis Tissue Forcep, 150mm, 4/5 teeth
- 18910-50 Allis Tissue Forcep, 150mm, 5/6 teeth
- 18910-60 Allis Tissue Forcep, 150mm, 6/7 teeth



-
- 18920-16 Babcock Tissue Forcep, serrated jaws, 165mm
 - 18920-20 Babcock Tissue Forcep, serrated jaws, 200mm
 - 18920-23 Babcock Tissue Forcep, serrated jaws, 230mm




-
- 18930-13 Duval Tissue Forcep, serrated, 190mm, 13mm jaws



-
- 18940-19 Judd-Allis Tissue Forcep, 195mm, 3/4 teeth




General Tissue Forceps



18950-12 Lanes Tissue Forcep, 1/2 teeth, fenestrated, 125mm


18950-15 Lanes Tissue Forcep, 1/2 teeth, fenestrated, 150mm

18950-19 Lanes Tissue Forcep, 1/2 teeth, fenestrated, 195mm



18960-18 Littlewood Tissue Forcep, 2/3 teeth, 190mm


Tongue Forceps



19000-13 Guys Tongue Forcep, screw joint, child, 130mm

19000-17 Guys Tongue Forcep, screw joint, adult, 175mm

Towel Clips




19100-09 Backhaus Towel Clip, 90mm




HEALTH MISSION
INSTRUMENTS

General Towel Clips


19100-12 Backhaus Towel Clip, 125mm



19130-09 Frenchay Towel Clip, cross action, non-penetrating, 90mm




19150-09 Shackle Towel Clip, cross action, 90mm



Trocars

19170-16 Ascites Trocar, 16FG


19170-21 Ascites Trocar, 21FG



19200-08 Hydrocele Trocar, 8FG

19200-09 Hydrocele Trocar, 9FG

19200-10 Hydrocele Trocar, 10FG



19210-22 Nelson Trocar, 22FG


19210-27 Nelson Trocar, 27FG

19210-33 Nelson Trocar, 33FG



HEALTH MISSION
INSTRUMENTS

General Vein Strippers




- 19280-06 Graves Myers Vein Stripper, acorn only, 6mm
- 19280-10 Graves Myers Vein Stripper, acorn only, 10mm
- 19280-13 Graves Myers Vein Stripper, acorn only, 13mm
- 19280-50 Graves Myers Vein Stripper, cable only, 900mm
- 19280-60 Graves Myers Vein Stripper, 'T' handle only
- 19280-90 Graves Myers Vein Stripper, complete, 3 acorns, cable, 'T' handle



- 19289-05 Mayo External Vein Stripper, 330mm, 5mm bore, 8mm o.d., Codman pattern
- 19289-07 Mayo External Vein Stripper, 330mm, 7mm bore, 10mm o.d., Codman pattern
- 19289-09 Mayo External Vein Stripper, 330mm, 9mm bore, 12mm o.d., Codman pattern



- 19290-05 Mayo External Vein Stripper, 385mm, 5mm bore, 8mm o.d., Pilling pattern
- 19290-07 Mayo External Vein Stripper, 385mm, 7mm bore, 10mm o.d., Pilling pattern
- 19290-09 Mayo External Vein Stripper, 385mm, 9mm bore, 12mm o.d., Pilling pattern




- 19300-06 Myer Vein Stripper, length 890mm, 6mm
- 19300-10 Myer Vein Stripper, length 890mm, 10mm
- 19300-13 Myer Vein Stripper, length 890mm, 13mm
- 19300-50 Myer Vein Stripper, handle only



General Vein Strippers

- 19310-06 Nabatoff Vein Stripper, acorn only, 6mm, 7BA
- 19310-09 Nabatoff Vein Stripper, acorn only, 9mm, 7BA
- 19310-12 Nabatoff Vein Stripper, acorn only, 12mm, 7BA
- 19310-15 Nabatoff Vein Stripper, acorn only, 15mm, 7BA
- 19310-20 Nabatoff Vein Stripper, cylindrical probe tip only, 10swg, 7BA
- 19310-50 Nabatoff Vein Stripper, cable only, 16 swg, 900mm, 7BA
- 19310-55 Nabatoff Vein Stripper, rack only, 7BA
- 19310-60 Nabatoff Vein Stripper, 'T' handle only, 7BA
- 19310-65 Nabatoff Vein Stripper, case only
- 19310-90 Nabatoff Vein Stripper, complete in case, 4 acorns, 2 cylindrical probe tips, cable, 'T' handle, rack, 7BA
- 19320-06 Nabatoff Vein Stripper, acorn only, 6mm, 6BA
- 19320-09 Nabatoff Vein Stripper, acorn only, 9mm, 6BA
- 19320-12 Nabatoff Vein Stripper, acorn only, 12mm, 6BA
- 19320-15 Nabatoff Vein Stripper, acorn only, 15mm, 6BA
- 19320-20 Nabatoff Vein Stripper, cylindrical probe tip only, 10swg, 6BA
- 19320-50 Nabatoff Vein Stripper, cable only, 16 swg, 900mm, 6BA
- 19320-55 Nabatoff Vein Stripper, rack only, 6BA
- 19320-60 Nabatoff Vein Stripper, 'T' handle only, 6BA



HEALTH MISSION
INSTRUMENTS




21000-22 Obwegeser Alveolar-Perinasal Wiring Awl, 225mm



21010-10 Kelsey Fry Bone Awl, curved, standard pattern

21015-10 Kelsey Fry Bone Awl, curved, lightweight pattern



21020-14 Kelsey Fry Bone Awl, straight, standard pattern, 145mm

21025-14 Kelsey Fry Bone Awl, straight, lightweight pattern, 145mm



21030-14 Obwegeser Mandibular Awl, 140mm



21040-18 Obwegeser Zygomatic Suspension Awl, 180mm

Chisels



21200-05 Braithwaite Chisel, guard on left, 5mm

21210-05 Braithwaite Chisel, guard on right, 5mm

21220-05 Braithwaite Chisel, unguarded, 5mm



21230-10 Cottle Nasal Chisel, No. 1, 5mm, 185mm



21230-20 Cottle Nasal Chisel, No. 2, 6mm, 185mm




21230-30 Cottle Nasal Chisel, No. 3, 7mm, 185mm



21230-40 Cottle Nasal Chisel, No. 4, 9mm, 185mm



21230-50 Cottle Nasal Chisel, No. 5, 11mm, 185mm



21230-60 Cottle Nasal Chisel, No. 6, 16mm, 185mm, fish tail



Plastic and Maxillo Facial Chisels

- 21235-03 Eastman Dental Chisel, length 190mm, 3mm
21235-05 Eastman Dental Chisel, length 190mm, 5mm



- 21238-04 Fomon Chisel, with guard, 165mm, 4mm
21238-05 Fomon Chisel, with guard, 165mm, 5mm
21238-06 Fomon Chisel, with guard, 165mm, 6mm
21238-07 Fomon Chisel, with guard, 165mm, 7mm



- 21240-18 Frenchay Chisel, curved to left, length 180mm



- 21250-18 Frenchay Chisel, curved to right, length 180mm




- 21260-18 Frenchay Chisel, straight, length 180mm



- 21270-02 Jenkins Chisel, length 145mm, 2mm
21270-04 Jenkins Chisel, length 145mm, 4mm
21270-06 Jenkins Chisel, length 145mm, 6mm
21270-08 Jenkins Chisel, length 145mm, 8mm
21270-10 Jenkins Chisel, length 145mm, 10mm
21270-12 Jenkins Chisel, length 145mm, 12mm
21270-14 Jenkins Chisel, length 145mm, 14mm



- 21280-10 MacEwen Chisel, hard edge, 10mm
21280-13 MacEwen Chisel, hard edge, 13mm
21280-20 MacEwen Chisel, hard edge, 20mm
21280-25 MacEwen Chisel, hard edge, 25mm
21280-30 MacEwen Chisel, hard edge, 30mm



- 21290-09 McIndoe Chisel, length 150mm, 9mm
21290-11 McIndoe Chisel, length 150mm, 11mm
21290-13 McIndoe Chisel, length 150mm, 13mm
21290-15 McIndoe Chisel, length 150mm, 15mm



Plastic and Maxillo Facial Chisels



21300-09 McIndoe Chisel, mushroom head, length 150mm, 9mm

21300-11 McIndoe Chisel, mushroom head, length 150mm, 11mm

21300-13 McIndoe Chisel, mushroom head, length 150mm, 13mm

21300-15 McIndoe Chisel, mushroom head, length 150mm, 15mm



21310-02 Read Chisel, length 180mm, 2mm

21310-03 Read Chisel, length 180mm, 3mm

21310-04 Read Chisel, length 180mm, 4mm

21310-05 Read Chisel, length 180mm, 5mm

21310-06 Read Chisel, length 180mm, 6mm

21310-07 Read Chisel, length 180mm, 7mm

21310-08 Read Chisel, length 180mm, 8mm




21320-03 Ward Chisel, length 180mm, 3mm

21320-04 Ward Chisel, length 180mm, 4mm

21320-05 Ward Chisel, length 180mm, 5mm

Collar Crushing Forceps



21400-14 Denis Browne Collar Crushing Forceps, 140mm

Curettes



21500-17 Exner Curette, double ended, angled oval ends, 175mm



21510-17 Exner Curette, double ended, angled round ends, 175mm



21520-17 Exner Curette, double ended, oval/round ends, 175mm



HEALTH MISSION
INSTRUMENTS

Plastic and Maxillo Facial Curettes

21530-18 McEwen Curette, double ended, 180mm, oval



21540-18 McEwen Curette, double ended, 180mm, round



Dissecting Forceps

21700-10 Mount Vernon Dissecting Forcep, serrated, 100mm



21705-14 Lanes Dissecting Forcep, rough handle, 1/2 teeth, 140mm




21705-18 Lanes Dissecting Forcep, rough handle, 1/2 teeth, 180mm




Elevators

21800-23 Bristow Malar Elevator, 230mm




21810-10 Cryer Exolever, left



21820-10 Cryer Exolever, right




Plastic and Maxillo Facial Elevators



21830-10 Cryer Exclever, straight



21840-22 Howarth Elevator, double ended, 220mm




21850-18 Kilner Malar Elevator, 180mm



21853-21 Malar Elevator, for elevating the zygomatic arch and malar bone by intra-oral approach, 210mm



21855-16 Mitchell Trimmer, double ended, 160mm



21860-21 Rowe Zygomatic Elevator



21863-17 Roberts Periosteal Elevator, 170mm



21865-15 Read Periosteal Elevator, single ended, blunt, 150mm



Plastic and Maxillo Facial Elevators

- 21866-21 Read Periosteal Elevator, double ended, blunt, 215mm



- 21870-17 Syme Periosteal Elevator, double ended, 170mm



- 21880-17 Ward Periosteal Elevator, fan shape, 170mm



- 21890-10 Warwick James Elevator, left



- 21900-10 Warwick James Elevator, right




- 21910-10 Warwick James Elevator, straight



Files

- 22000-17 Alveolectomy File, double ended, 170mm



Gouges

- 22170-03 Coupland Gouge, 3mm
22170-04 Coupland Gouge, 4mm
22170-05 Coupland Gouge, 5mm



- 22195-10 Parkes Hump Gouge, double guarded



Plastic and Maxillo Facial Gouges



- 22200-02 Read Gouge, length 180mm, 2mm
- 22200-03 Read Gouge, length 180mm, 3mm
- 22200-04 Read Gouge, length 180mm, 4mm
- 22200-05 Read Gouge, length 180mm, 5mm
- 22200-06 Read Gouge, length 180mm, 6mm

Hand Drills



- 22300-19 Archimedean Hand Drill, 3mm fitting, 190mm

Hooks



- 22380-15 Harrison Tendon Hook, 150mm



- 22400-21 Poswillo Malar Hook, 215mm


Mallets



- 22500-06 Fickling Dental Mallet, s/s, 190mm, 25mm diameter, 190g (6.3/4oz)




- 22505-11 Hammer For Hand Surgery, s/s, 195mm, 320g (11.1/4oz)



- 22508-17 Mallet For Ultra Fine Osteotome, s/s, nylon face insert, tuftol handle, 170mm, 20mm diameter, 145g (5oz)



HEALTH MISSION
INSTRUMENTS



22800-12 Joseph Nasal Saw, horizontal shaft, 125mm x 85mm, cranked to the left

22810-12 Joseph Nasal Saw, horizontal shaft, 125mm x 85mm, cranked to the right




22820-18 Joseph Nasal Saw, straight, 180mm




22830-18 Joseph Nasal Saw, vertical shaft, 185mm, cranked to the left

22840-18 Joseph Nasal Saw, vertical shaft, 185mm, cranked to the right



22860-20 Mowlem Nasal Saw, 200mm

Needleholders




22900-12 Foster Gillies Needleholder, left hand, 130mm

22910-12 Foster Gillies Needleholder, right hand, 130mm



22920-16 Gillies Needleholder, left hand, 165mm


22930-16 Gillies Needleholder, right hand, 165mm



Needleholders - Tungsten Carbide


23000-13 Foster Gillies Needleholder, tungsten carbide, left hand, 130mm

23010-13 Foster Gillies Needleholder, tungsten carbide, right hand, 130mm



23020-16 Gillies Needleholder, tungsten carbide, left hand, 165mm


23030-16 Gillies Needleholder, tungsten carbide, right hand, 165mm



Oral Forceps


23090-18 Fickling Oral Forcep, 1/2 teeth, 180mm

23091-18 Fickling Oral Forcep, serrated, 180mm





23110-23 Rowe Disimpaction Oral Forcep, left, screw joint, 230mm



23120-23 Rowe Disimpaction Oral Forcep, right, screw joint, 230mm

Osteotomes



23180-10 Cinelli Osteotome, double guarded, 10mm
23180-12 Cinelli Osteotome, double guarded, 12mm
23180-14 Cinelli Osteotome, double guarded, 14mm
23180-16 Cinelli Osteotome, double guarded, 16mm




23200-02 French (Gillies) Osteotome, length 140mm, 2mm
23200-03 French (Gillies) Osteotome, length 140mm, 3mm
23200-04 French (Gillies) Osteotome, length 140mm, 4mm
23200-05 French (Gillies) Osteotome, length 140mm, 5mm
23200-06 French (Gillies) Osteotome, length 140mm, 6mm
23200-07 French (Gillies) Osteotome, length 140mm, 7mm
23200-08 French (Gillies) Osteotome, length 140mm, 8mm
23200-09 French (Gillies) Osteotome, length 140mm, 9mm
23200-10 French (Gillies) Osteotome, length 140mm, 10mm
23200-11 French (Gillies) Osteotome, length 140mm, 11mm
23200-12 French (Gillies) Osteotome, length 140mm, 12mm
23200-13 French (Gillies) Osteotome, length 140mm, 13mm



23210-02 Jayes Osteotome, 2mm
23210-04 Jayes Osteotome, 4mm
23210-05 Jayes Osteotome, 5mm
23210-07 Jayes Osteotome, 7mm
23210-10 Jayes Osteotome, 10mm
23210-12 Jayes Osteotome, 12mm
23210-15 Jayes Osteotome, 15mm




- 23215-11 MacGregor Nasal Osteotome, bevelled corners, flat handle with mushroom head, 11mm
- 23215-13 MacGregor Nasal Osteotome, bevelled corners, flat handle with mushroom head, 13mm
- 23215-15 MacGregor Nasal Osteotome, bevelled corners, flat handle with mushroom head, 15mm




- 23220-11 Matthew Osteotome, bevelled at one corner, 11mm



- 23226-10 Neivert Osteotome, single guarded, left, 3.5mm
- 23226-20 Neivert Osteotome, single guarded, right, 3.5mm
- 23226-30 Neivert Osteotome, single guarded, straight, 3.5mm



- 23230-02 Read Osteotome, length 180mm, 2mm
- 23230-03 Read Osteotome, length 180mm, 3mm
- 23230-04 Read Osteotome, length 180mm, 4mm
- 23230-05 Read Osteotome, length 180mm, 5mm
- 23230-06 Read Osteotome, length 180mm, 6mm



- 23240-10 Rozner Osteotome, probe point on each side of inclined cutting edge





- 23243-10 Ultra Fine Osteotome, case only
- 23243-20 Ultra Fine Osteotome, rack only, to take 1 x mallet and 6 x osteotomes
- 23243-90 Ultra Fine Osteotome, case and rack, complete, to take 1 x mallet and 6 x osteotomes
- 23244-02 Ultra Fine Osteotome, TiN coated blade, tufnol handle, 140mm, 2mm
- 23244-04 Ultra Fine Osteotome, TiN coated blade, tufnol handle, 140mm, 4mm
- 23244-06 Ultra Fine Osteotome, TiN coated blade, tufnol handle, 140mm, 6mm
- 23244-08 Ultra Fine Osteotome, TiN coated blade, tufnol handle, 140mm, 8mm
- 23244-10 Ultra Fine Osteotome, TiN coated blade, tufnol handle, 140mm, 10mm
- 23244-12 Ultra Fine Osteotome, TiN coated blade, tufnol handle, 140mm, 12mm
- 23244-90 Ultra Fine Osteotome, TiN coated blade, tufnol handle, 140mm, set of 6 (2mm to 12mm)
- 23245-02 Ultra Fine Osteotome, TiN coated blade, tufnol handle, 160mm, 2mm
- 23245-04 Ultra Fine Osteotome, TiN coated blade, tufnol handle, 160mm, 4mm
- 23245-06 Ultra Fine Osteotome, TiN coated blade, tufnol handle, 160mm, 6mm
- 23245-08 Ultra Fine Osteotome, TiN coated blade, tufnol handle, 160mm, 8mm
- 23245-10 Ultra Fine Osteotome, TiN coated blade, tufnol handle, 160mm, 10mm
- 23245-12 Ultra Fine Osteotome, TiN coated blade, tufnol handle, 160mm, 12mm
- 23245-90 Ultra Fine Osteotome, TiN coated blade, tufnol handle, 160mm, set of 6 (2mm to 12mm)



- 23250-03 Ward Osteotome, length 180mm, 3mm
- 23250-04 Ward Osteotome, length 180mm, 4mm
- 23250-05 Ward Osteotome, length 180mm, 5mm

Raspatories



- 23300-10 Cleft Palate Raspatory, curved downwards
- 23310-10 Cleft Palate Raspatory, curved upwards



- 23320-10 Cleft Palate Raspatory, curved upwards, left



- 23330-10 Cleft Palate Raspatory, curved upwards, right




- 23350-10 Joseph Nasal Raspatory, curved
23360-10 Joseph Nasal Raspatory, straight




Rasps

- 23400-10 Barsky Glabella Rasp, curved, reverse cutting teeth, concave surface



- 23410-10 Barsky Nasal Rasp, straight, reverse cutting teeth, concave surface



- 23435-17 Maltz Nasal Rasp, straight, with longitudinal serrations, 170mm



- 23440-19 Maltz Nasal Rasp, curved, deep cut serrations, 190mm



- 23450-17 Maltz Nasal Rasp, straight, deep cut serrations, 175mm

- 23455-17 McIndoe Nasal Rasp, straight, reverse cutting teeth, 170mm




- 23460-19 Reidy Nasal Rasp, curved, cross serrated, 195mm






- 23490-02 Harrison Nicholl Reamer, s/s, double ended, malleable, size 1/2
- 23490-04 Harrison Nicholl Reamer, s/s, double ended, malleable, size 3/4
- 23490-06 Harrison Nicholl Reamer, s/s, double ended, malleable, size 5/6
- 23490-08 Harrison Nicholl Reamer, s/s, double ended, malleable, size 7/8
- 23490-10 Harrison Nicholl Reamer, s/s, double ended, malleable, size 9/10
- 23490-12 Harrison Nicholl Reamer, s/s, double ended, malleable, size 11/12
- 23490-14 Harrison Nicholl Reamer, s/s, double ended, malleable, size 13/14
- 23490-90 Harrison Nicholl Reamer, s/s, double ended, malleable, set of 7 sizes (1/2 - 13/14)


Retractors



- 23500-12 Awty Osteotomy Retractor, 12mm wide, pair
- 23500-17 Awty Osteotomy Retractor, 17mm wide, pair




- 23505-25 Attwood Retractor, with through handle fibrelight and suction tube, 2mm thick blade, 110mm deep x 25mm wide
- 23505-40 Attwood Retractor, with through handle fibrelight and suction tube, 2mm thick blade, 150mm deep x 40mm wide
- 23505-50 Attwood Retractor, with through handle fibrelight and suction tube, 2mm thick blade, 190mm deep x 50mm wide



- 23510-19 Bowdler-Henry (Monks Tissue) Rake Retractor, 170mm



- 23515-25 Converse Retractor, with fibrelight, 25mm deep x 12mm wide



- 23520-10 Cottle Upper Lateral Exposing Retractor



- 23530-08 Desmarres Retractor, tight curve, 140mm, 8mm wide

- 23530-12 Desmarres Retractor, tight curve, 140mm, 12mm wide


- 23530-14 Desmarres Retractor, tight curve, 140mm, 14mm wide

- 23530-16 Desmarres Retractor, tight curve, 140mm, 16mm wide


- 23530-18 Desmarres Retractor, tight curve, 140mm, 18mm wide




- 23532-30 Eddey Parotid Retractor, dished blade, serrated, 30mm deep x 20mm wide blade



- 23534-17 Fitzpatrick Mandibular 'Trowel' Retractor, 170mm







- 23536-10 Fitzpatrick Tongue Guard Retractor, malleable, No.1, 84mm, suitable for small mouth or child
- 23536-20 Fitzpatrick Tongue Guard Retractor, malleable, No.2, 88mm, universal, does not need reshaping for opposite sides of the mouth
- 23536-30 Fitzpatrick Tongue Guard Retractor, malleable, No.3, 86mm, a variation in contour for No.2
- 23536-40 Fitzpatrick Tongue Guard Retractor, malleable, No.4, 88mm, suitable for large or edentulous mouth
- 23536-90 Fitzpatrick Tongue Guard Retractor, malleable, set of 4



- 23538-18 Freeman Face Lift Retractor, 185mm



- 23540-16 Hajeck Retractor, for upper lip, 165mm



- 23545-15 Henderson Stripping Retractor, for ramus, 150mm




Plastic and Maxillo Facial Retractors

- 23547-16 Intra Oral Vertical Osteotomy Retractor, with through handle fibrelight, 100mm deep x 16mm wide




- 23550-14 Kay Austin Retractor, 140mm



- 23560-15 Kilner 'Cats Paw' Retractor, double ended, 150mm



- 23565-25 Hurren Retractor, with through handle fibrelight, 60 deg blade angle, without lip, 150mm deep x 25mm wide




- 23570-10 Kilner Alae Retractor, blunt, 85mm x 10mm
23570-13 Kilner Alae Retractor, blunt, 85mm x 13mm
23570-15 Kilner Alae Retractor, blunt, 85mm x 15mm
23580-10 Kilner Alae Retractor, sharp, 85mm x 10mm
23580-13 Kilner Alae Retractor, sharp, 85mm x 13mm
23580-15 Kilner Alae Retractor, sharp, 85mm x 15mm







- 23590-18 Kilner Cheek Retractor, double ended, 25mm x 18mm
- 23590-22 Kilner Cheek Retractor, double ended, 25mm x 22mm
- 23590-27 Kilner Cheek Retractor, double ended, 35mm x 27mm
- 23590-28 Kilner Cheek Retractor, double ended, 30mm x 28mm
- 23600-18 Kilner Cheek Retractor, double ended, insulated, 25mm x 18mm
- 23600-22 Kilner Cheek Retractor, double ended, insulated, 25mm x 22mm
- 23600-27 Kilner Cheek Retractor, double ended, insulated, 35mm x 27mm
- 23600-28 Kilner Cheek Retractor, double ended, insulated, 30mm x 28mm



- 23610-15 Kilner Nasal Retractor, 150mm




- 23611-05 Kilner Skin Retractor, long teeth, 50mm




- 23611-06 Kilner Skin Retractor, small teeth, 60mm




- 23612-25 Mammary Retractor, with through handle fibrelight and suction tube, plain lip, 90mm deep x 25mm wide
- 23615-25 Mammary Retractor, with through handle fibrelight and suction tube, toothed lip, 150mm deep x 25mm wide, 5 teeth, semi-sharp
- 23615-32 Mammary Retractor, with through handle fibrelight and suction tube, toothed lip, 150mm deep x 32mm wide, 6 teeth, semi-sharp
- 23615-40 Mammary Retractor, with through handle fibrelight and suction tube, toothed lip, 150mm deep x 40mm wide, 8 teeth, semi-sharp
- 23616-25 Mammary Retractor, with through handle fibrelight and suction tube, plain lip, 150mm deep x 25mm wide
- 23616-40 Mammary Retractor, with through handle fibrelight and suction tube, plain lip, 150mm deep x 40mm wide
- 23617-25 Mammary Retractor, with through handle fibrelight and suction tube, toothed lip, 190mm deep x 25mm wide, 5 teeth, semi-sharp
- 23618-25 Mammary Retractor, with through handle fibrelight and suction tube, plain lip, 190mm deep x 25mm wide
- 23623-25 Mammary Retractor, with through handle fibrelight and suction tube, plain lip, 110mm deep x 25mm wide
- 23623-40 Mammary Retractor, with through handle fibrelight and suction tube, plain lip, 110mm deep x 40mm wide




-
- 23634-25 Mammary Retractor, with through handle fibrelight, without lip, 150mm deep x 25mm wide
 - 23639-40 Mammary Retractor, with through handle fibrelight, without lip, 190mm deep x 40mm wide
 - 23639-50 Mammary Retractor, with through handle fibrelight, without lip, 190mm deep x 50mm wide
 - 23644-25 Mammary Retractor, with through handle fibrelight, without lip, 110mm deep x 25mm wide





- 23631-25 Mammary Retractor, with through handle fibrelight, toothed lip, 90mm deep x 25mm wide, 5 teeth, semi-sharp
- 23631-40 Mammary Retractor, with through handle fibrelight, toothed lip, 90mm deep x 40mm wide, 8 teeth, semi-sharp
- 23635-25 Mammary Retractor, with through handle fibrelight, toothed lip, 150mm deep x 25mm wide, 5 teeth, semi-sharp
- 23635-40 Mammary Retractor, with through handle fibrelight, toothed lip, 150mm deep x 40mm wide, 8 teeth, semi-sharp
- 23637-25 Mammary Retractor, with through handle fibrelight, toothed lip, 190mm deep x 25mm wide, 5 teeth, semi-sharp
- 23637-40 Mammary Retractor, with through handle fibrelight, toothed lip, 190mm deep x 40mm wide, 8 teeth, semi-sharp
- 23642-25 Mammary Retractor, with through handle fibrelight, toothed lip, 110mm deep x 25mm wide, 5 teeth, semi-sharp
- 23642-40 Mammary Retractor, with through handle fibrelight, toothed lip, 110mm deep x 40mm wide, 8 teeth, semi-sharp
- 23642-44 Mammary Retractor, with through handle fibrelight, toothed lip, 110mm deep x 44mm wide, 9 teeth, semi-sharp



- 23632-25 Mammary Retractor, with through handle fibrelight, plain lip, 90mm deep x 25mm wide
- 23632-40 Mammary Retractor, with through handle fibrelight, plain lip, 90mm deep x 40mm wide
- 23636-25 Mammary Retractor, with through handle fibrelight, plain lip, 150mm deep x 25mm wide
- 23636-40 Mammary Retractor, with through handle fibrelight, plain lip, 150mm deep x 40mm wide
- 23638-25 Mammary Retractor, with through handle fibrelight, plain lip, 190mm deep x 25mm wide
- 23638-40 Mammary Retractor, with through handle fibrelight, plain lip, 190mm deep x 40mm wide
- 23643-25 Mammary Retractor, with through handle fibrelight, plain lip, 110mm deep x 25mm wide
- 23643-40 Mammary Retractor, with through handle fibrelight, plain lip, 110mm deep x 40mm wide
- 23643-44 Mammary Retractor, with through handle fibrelight, plain lip, 110mm deep x 44mm wide



- 23640-10 Laster Cheek Retractor, for third upper molar, enables a chisel to be held, 100mm



Plastic and Maxillo Facial Retractors


- 23648-13 McIndoe Retractor, temporo-mandibular, 13mm wide x 25mm deep




- 23650-44 Oringer Cheek And Lip Retractor, size 1, 44mm deep
23650-51 Oringer Cheek And Lip Retractor, size 2, 51mm deep
23650-57 Oringer Cheek And Lip Retractor, size 3, 57mm deep




- 23655-10 Mustarde Offset Orbital Retractor, central, 110mm




- 23655-20 Mustarde Offset Orbital Retractor, deep right, 110mm






23655-30 Mustarde Offset Orbital Retractor, shallow right, 110mm



23655-40 Mustarde Offset Orbital Retractor, deep left, 110mm



23655-50 Mustarde Offset Orbital Retractor, shallow left, 110mm



23660-17 Prognathism Channel Retractor, 170mm

